

BUTTERFLY IQ3 ULTRASOUND ENTERPRISE BUNDLE

Learn more about this and other Products with a Mission® at www.NARescue.com



UNPARALLELED IMAGING QUALITY WITH LIGHTNING-FAST DATA PROCESSING

Meet the new Butterfly iQ3™ - our easiest-to-use point-of-care ultrasound (POCUS) device with the sharpest image quality, rapid data processing, AI, advanced imaging tools, user-centric ergonomics, and outstanding scan time. In the field, you often have only minutes to make difficult diagnostic decisions that can impact patients' survival. What if you could have the power of whole-body ultrasound with you?

The Butterfly iQ3™ is the clear answer in point-of-care ultrasound. The iQ3 has superior data transfer rate, 2x faster processing power, user-centric ergonomics and configurable buttons. The user-friendly controls allow for seamless navigation, enabling quick scans and interpretations, thereby saving critical time during diagnoses. Designed to make ultrasound easy to use. Simply connect the probe to your compatible iOS or Android device, select from 20+ presets in seconds with the swipe of a finger — and start scanning.

FEATURES/SPECIFICATIONS

- The clear answer in point-of-care ultrasound
- Operating software compatible with iOS and Android
- Digital image storage of 1 minute. Unlimited digital cloud storage
- 2220 mAh Lithium Ion battery (Not replaceable)
- Contact charger included. Battery Life Up to 2 hours
- Battery recharge from 10% to 90% in 90 minutes
- Single probe solution covering 20+ imaging applications
- Ultra portable device, can be used inside and outside the rig
- For prescription use only. Designed for external ultrasound imaging
- Durable, military-grade tested hardware (MIL-STD-810G)
- 3 year manufacturer's warranty on device
- Bundle includes ScanLab, Butterfly Academy, and On-Demand Videos for lifetime of device (no recurring fees with bundle option)
- Purchasers will receive access to ScanLab training system (ScanLab is iPad compatible only)

IQ3™ SPECIFICATIONS

- **Frequency:** 1-12 MHz
- **Butterfly ultrasound chip:** Next generation
- **Data transfer rate:** 9.6Gb/s
- **Frame rate:**
 - » 28Hz (Cardiac*1)
 - » 28Hz (Abdomen)
- **Probe Body:** 152 mm x 52 mm x 37 mm
- **Weight distribution:** Ergonomically distributed (300g)
- **Probe dimensions:**
 - » Area: 808 mm²
 - » Length: 35.3 mm
 - » Width: 22.9 mm
- **Buttons:** 3 buttons; configurable
- **Charging:** Magnetic 2-hour contact charging
- **Standard warranty:** 3 years



DIMENSIONS

- Packaged: H 6.4 in. x W 2.2 in. x D 1.4 in.
- Weight: 0.68 lbs.

20+ CLINICAL PRESETS

- Abdomen
- Abdomen Deep
- Aorta & Gallbladder
- Bladder
- Cardiac (Cardiac Deep, Coherence Imaging)
- FAST
- Lung (Lung Tissue)
- MSK-Soft Tissue
- Musculoskeletal
- Nerve
- OB 1/GYN
- OB 2/3
- Ophthalmic
- Pediatric Abdomen
- Pediatric Cardiac
- Pediatric Lung
- Small Organ
- Vascular: Access
- Vascular: Carotid
- Vascular: Deep Vein
- Vascular: Face

IMAGING MODES

- B-Mode
- M-Mode
- Color Doppler
- Power Doppler
- Pulsed-Wave Doppler
- Biplane Imaging™

BUTTERFLY ULTRASOUND

ITEM #	DESCRIPTION
20-0123	Butterfly iQ3 Ultrasound Enterprise Bundle
ZZ-1829	Butterfly Ultrasound - Hard Case Only
ZZ-1830	Butterfly Ultrasound - Soft Case Only
ZZ-1831	Butterfly Ultrasound - Soft Case Mini Only



* Phone not included.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277



BUTTERFLY ULTRASOUND HARD CASE | ZZ-1829

Under a durable shell, the redesigned Butterfly Hard Case secures your probe, charger, full bottle of gel, and mobile device. Neatly loop your cord around a custom cable wrap.

Secure your case with a satisfying click, and lock it away for travel with two eyelets. Take it on the go with a new ergonomic handle and a lid which doubles as a mobile device stand.

DIMENSIONS

- H 10.8 in. x W 10.4 in. x D 3.5 in.
- Weight: 3.4 lb.

MATERIALS

- Plastic: Poly carbonate/acrylonitrile butadiene styrene (PC-ABS)
- Handle: Thermoplastic (TPU)
- Inner Tray: Silicone

SPECIFICATIONS

- Durable shell compatible with the iQ+ and iQ3™ Butterfly Ultrasounds
- Ergonomic handle, with a lid that can double as a mobile device stand
- The hard case is made of poly carbonate/acrylonitrile butadiene styrene (PC-ABS), while the inner tray is made of silicone.
- Case only. Does not include probe, charger, tablet, or gel



BUTTERFLY ULTRASOUND SOFT CASE | ZZ-1830

The Butterfly Soft Case is lightweight, durable and makes it easy to store and carry your Butterfly probe, charger, and mobile device. It even has a shoulder strap and room for your gel.

DIMENSIONS

- H 8.3 in. x W 7 in. x D 1.5 in.
- Weight: 0.62 lb.

SPECIFICATIONS

- Lightweight, durable case compatible with the iQ+ and iQ3™ Butterfly Ultrasounds
- Case is made of 1050D Ballistic Nylon
- Case only. Does not include probe, charger, mobile device, or gel



* Case only. Device not included.

* Case only. Device not included.



BUTTERFLY ULTRASOUND SOFT CASE MINI | ZZ-1831

The Butterfly Soft Case Mini is a lightweight and compact solution for carrying your probe, charger, and a small tube of gel wherever you go. This sleek, portable case effortlessly travels with you, ensuring your probe, charger, and a small tube of gel are always at hand.

DIMENSIONS

- H 8.9 in. x W 3.1 in. x D 1.9 in.
- Weight: 0.3 lb.

SPECIFICATIONS

- Lightweight, durable case compatible with the iQ+ and iQ3™ Butterfly Ultrasounds
- Case is made of 600D Polyester
- Case only. Does not include probe, charger, or gel



* Case only. Device not included.

BUTTERFLY ULTRASOUND

ITEM #	DESCRIPTION
20-0123	Butterfly iQ3 Ultrasound Enterprise Bundle
ZZ-1829	Butterfly Ultrasound - Hard Case Only
ZZ-1830	Butterfly Ultrasound - Soft Case Only
ZZ-1831	Butterfly Ultrasound - Soft Case Mini Only



Learn more about this and other Products with a Mission® at www.NARescue.com