

TRAUMA AND FIRST AID KIT WALL MOUNT HOOK

Learn more about this and
other *Products with a Mission*[®]
at www.NARescue.com



Designed To Hold North American Rescue's Trauma And First Aid Kits

North American Rescue's Trauma and First Aid Kit Wall Mount Hook provides high visibility and rapid access to critical medical supplies for use by initial first responders and employees in worksite/construction areas in the workplace.

The heavy-duty zinc plated steel hook is designed to hold North American Rescue's Trauma and First Aid Kits by their handles. The hook mounts to any wall, features rounded corners help prevent injury if bumped into, and comes complete with mounting hardware. To best protect your critical lifesaving products, hooks should be mounted indoors out of the direct elements.



FEATURES/SPECIFICATIONS

- Designed to store Class A or B Trauma and First Aid Kits (sold separately) in areas where rapid access is needed
- Heavy duty zinc plated steel hook
- Easily mounts to any wall

TFAK WALL MOUNT

ITEM #

ZZ-1819

DIMENSIONS

- 5/16 in. thick Steel
- Hook Dimensions - H 2.87 in. x W 3 in. x D 7.5 in
- Mounting Plate Dimensions: H 3 in. x W 5 in. x D 0.12 in.
- Weight: 1 lb.



NORTH AMERICAN RESCUE[®]
www.NARescue.com • 888.689.6277