

PUBLIC ACCESS INDIVIDUAL BLEEDING CONTROL KIT

VACUUM-SEALED

Control the bleed with these easy-to-use tools!

Life threatening bleeding can occur from everyday emergencies such as vehicle crashes, incidents involving farming equipment or public transportation, accidental injuries from natural disasters like tornadoes, and inflicted injuries from man-made disasters such as shootings and bombings. The Individual Bleeding Control Kit by North American Rescue is designed to provide essential equipment that empowers the general public to take action as immediate responders in stopping life threatening bleeding.

BASIC:

The NAR Individual Bleeding Control Kit provides life saving bleeding control equipment such as tourniquets, pressure dressings and gauze bandages and puts them into the hands of the public to help control bleeding and save lives.

INTERMEDIATE:

Includes all equipment found in the Basic Individual Bleeding Control Kit with the added capability to treat penetrating chest wounds with the HyFin® Vent Compact Chest Seal Twin Pack (for entrance and exit wounds).

ADVANCED / ADVANCED BCD / ADVANCED CELOX™:

Includes all equipment found in the Intermediate Individual Bleeding Control Kit with the added capability to treat major bleeding with a hemostatic dressing to enhance blood clotting abilities for faster bleeding control. Contained in a ruggedized, vacuum-sealed pouch, this life saving kit provides intuitive and easy-to-use tools that are proven to help to save lives such as the fastest, safest and most effective prehospital tourniquet in the world, the C-A-T®, quality gauze, shears, and pressure dressings/bandages. Included are step-by-step illustrated "just in time" pictorial instructions.

FEATURES/SPECIFICATIONS

- Designed to provide essential equipment that empowers the general public to take action as immediate responders in stopping life threatening bleeding
- Provides intuitive and easy-to-use tools that are proven to help to save lives
- Kits are available in Basic, Intermediate, Advanced, Advanced BCD, and Advanced Celox™ configurations
- Packaged in a compact, ruggedized vacuum-sealed pouch



Learn more about this and other Products with a Mission® at www.NARescue.com

PABC IND. KIT - VAC

ITEM #	OPTION
80-0465	Basic
80-0466	Intermediate
80-0467	Advanced
80-0522	Advanced BCD
80-1569	Advanced with Celox™ Gauze

* Explore our full line of Public Access products at www.NAResceu.com

DIMENSIONS

- H 7 in. x W 5 in. x D 2.75 in.
- Weight: 14 -17.4 oz



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PUBLIC ACCESS INDIVIDUAL BLEEDING CONTROL KIT

VACUUM-SEALED

all KITS INCLUDE:



KIT CONTENTS

- 1x C-A-T® Tourniquet - Orange
- 1x 6 in. NAR® Responder ETD® (Emergency Trauma Dressing)
- 1x NAR® Wound Packing Gauze
- 1x NAR® Survival Blanket
- 2 pr Responder Nitrile Gloves, Lg
- 1x Responder Trauma Shears, 7.25 in.
- 1x Permanent Marker, Small
- 1x Just in Time Instruction Card

- +plus, by kit type**
- 1x NAR® Wound Packing Gauze
 - or -
 - 1x NAR® Wound Packing Gauze
 - 1x HyFin® Compact, Vent Twin Pack
 - or -
 - 1x Combat Gauze® LE, Z-Fold Hemostatic Bandage
 - 1x HyFin® Compact, Vent Twin Pack
 - or -
 - 1x QuikClot® Bleeding Control Dressing
 - 1x HyFin® Compact, Vent Twin Pack
 - or -
 - 1x 5 ft. Celox™ Hemostatic Trainer
 - 1x HyFin® Compact, Vent Twin Pack

+plus, by kit type:



PABC IND. KIT - VAC

ITEM #	OPTION
80-0465	Basic
80-0466	Intermediate
80-0467	Advanced
80-0522	Advanced BCD
80-1569	Advanced with Celox™ Gauze

* Explore our full line of Public Access products at www.NAResceu.com

Learn more about this and other Products with a Mission® at www.NAResceu.com

