

PUBLIC ACCESS INDIVIDUAL BLEEDING CONTROL KIT

NYLON

Train To Stop Life Threatening Bleeding And Save A Life!

Life threatening bleeding can occur from everyday emergencies such as vehicle crashes, incidents involving farming equipment or public transportation, accidental injuries from natural disasters like tornadoes, and inflicted injuries from man-made disasters such as shootings and bombings. The Individual Bleeding Control Kit by North American Rescue is designed to provide essential equipment that empowers the general public to take action as immediate responders in stopping life threatening bleeding.

Contained in a nylon carry bag, this life saving kit provides intuitive and easy-to-use tools that are proven to help to save lives such as the fastest, safest and most effective prehospital tourniquet in the world, the C-A-T®, quality gauze, shears, and pressure dressings/bandages. Included are step-by-step illustrated "just in time" pictorial instructions.

Click here to purchase the Public Access Individual Bleeding Control Kit - BAG ONLY

NAR Individual Bleeding Control Kits are available in Basic, Intermediate, Advanced, Advanced BCD, and Advanced Celox™ configurations.

BASIC:

The NAR Individual Bleeding Control Kit provides life saving bleeding control equipment such as tourniquets, pressure dressings and gauze bandages and puts them into the hands of the public to help control bleeding and save lives.

INTERMEDIATE:

Includes all equipment found in the Basic Individual Bleeding Control Kit with the added capability to treat penetrating chest wounds with the HyFin® Vent Compact Chest Seal Twin Pack (for entrance and exit wounds).

ADVANCED / ADVANCED BCD / ADVANCED CELOX™ :

Includes all equipment found in the Intermediate Individual Bleeding Control Kit with the added capability to treat major bleeding with a hemostatic dressing that enhances blood clotting abilities—for faster bleeding control.

FEATURES/SPECIFICATIONS

- Designed to provide essential equipment that empowers the general public to take action as immediate responders in stopping life threatening bleeding
- Provides intuitive and easy-to-use tools that are proven to help to save lives
- Kits are available in Basic, Intermediate, Advanced, Advanced BCD, and Advanced Celox™ configurations
- Packaged in a compact, red nylon carry bag

Learn more about this and other Products with a Mission® at www.NARescue.com



PABC IND KIT - NYLON

ITEM #	OPTION
80-0452	Basic
80-0453	Intermediate
80-0454	Advanced
80-0523	Advanced BCD
80-1570	Advanced with Celox™

* Explore our full line of Public Access products at www.NARescue.com

DIMENSIONS

- H 7 in. x W 5 in. x D 2.5 in.
- Weight: 15 - 174 oz

PUBLIC ACCESS INDIVIDUAL BLEEDING CONTROL KIT

NYLON

all KITS INCLUDE:




+plus, by kit type:



KIT CONTENTS

- 1x C-A-T® Tourniquet
- 1x 6 in. NAR® Responder ETD® (Emergency Trauma Dressing)
- 2x NAR® Wound Packing Gauze - or -
- 2x NAR® Wound Packing Gauze
- 1x HyFin® Compact, Vent Twin Pack - or -
- 1x Combat Gauze® LE, Z-Fold Hemostatic Bandage
- 1x HyFin® Compact, Vent Twin Pack - or -
- 1x QuikClot® Bleeding Control Dressing
- 1x HyFin® Compact, Vent Twin Pack - or -
- 1x 5 ft. Celox™ First Aid Z-Fold Hemostatic Dressing
- 1x HyFin® Compact, Vent Twin Pack
- 2 pr Responder Nitrile Gloves, Lg
- 1x Responder Trauma Shears, 7.25 in.
- 1x NAR® Emergency Survival Blanket - 56 in. x 96 in.
- 1x Permanent Marker, Small
- 1x Just in Time Instruction Card
- 1x PVC Bleeding Control Patch
- 1x Red Nylon Bag

+plus, by kit type



PABC IND. KIT - NYLON

ITEM #	OPTION
80-0452	Basic
80-0453	Intermediate
80-0454	Advanced
80-0523	Advanced BCD
80-1570	Advanced with Celox™

* Explore our full line of Public Access products at www.NARRescue.com

Learn more about this and other Products with a Mission® at www.NARRescue.com

