

I-ROK XL

INDIVIDUAL RUGGED OUTDOOR KIT

Learn more about this and other **Products with a Mission®** at www.NARescue.com



Every good adventure requires that everyone gets home safely.

The I-ROK XL (Individual Rugged Outdoor Kit) from North American Rescue is the perfect companion for any outdoor environment. The I-ROK XL contains the same essential first-aid supplies of the Standard I-ROK Kit with extra supplies to accommodate multiple injuries. These kits are designed to treat common injuries including basic injuries such as cuts, scrapes, bites & stings as well as the life-saving equipment to control major bleeding.

The I-ROK XL comes in the larger ultimate personal R60 utility case that is crushproof, dustproof and adventure proof! Your first-aid and trauma supplies will be protected from the outdoor elements to help ensure your kit is ready if ever needed. Waterproof up to a depth of 2 meters for up to 30 minutes, IP68 rated protection from water, dirt, snow & dust and rubberized protective bumpers to protect against drops. The case handles can be used to secure the kit for your outing as you can strap it to your gear, ATV/UTV, motorcycle, camper or whatever you use to get you to where the adventures are. This ample hard case comes in a variety of colors and keeps your contents protected and secure.



Case handles allow you to open the kit while it's strapped to your gear, camper, motorcycle, ATV/UTV, etc.

FEATURES/SPECIFICATIONS

- Contains essential first-aid supplies to treat basic injuries as well as control major bleeding with extra supplies for multiple injuries
- Features a larger R60 storage case that is crushproof, dustproof and waterproof up to a depth of 2 meters
- Flexible lid organizer helps keep smaller medical supplies neatly stored
- Built-in pressure valve keeps water and dust out while balancing air pressure
- Larger hard case comes in a variety of colors and with the option of Basic or with Bleeding Control Dressing (hemostatic bandage) versions

I-ROK XL

ITEM #	COLOR	OPTION
82-0049	Black	Basic
82-0051	OD Green	Basic
82-0053	Black	w/ BCD
82-0055	OD Green	w/ BCD
82-0070	Black	w/ CELOX™
82-0072	OD Green	w/ CELOX™

DIMENSIONS

- Measurements (Interior): H 3 in. x W 8.9 in. x D 5.4 in.
- Measurements (Exterior): H 3.9 in. x W 11.1 in. x D 6.9 in.
- Weight: 3 lb



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277


I-ROK XL

INDIVIDUAL RUGGED OUTDOOR KIT

all KITS INCLUDE:



+plus, by kit type:



Learn more about this and other Products with a Mission® at www.NARescue.com



KIT CONTENTS

- 1x C-A-T® Tourniquet, Orange
- 1x Responder Compressed Gauze
 - or-
 - QuikClot® Hemostatic BCD
 - or-
 - CELOX™ First Aid Hemostatic Dressing
- 1x 4 in. NAR® Responder Mini ETD® (Emergency Trauma Dressing)
- 1x HyFin® Vent Chest Seal Twin Pack
- 10 x Butterfly Bandage
- 10 x Bandage (1 in. x 3 in.)
- 2 x Gauze Pad 2-pack (4 in. x 4 in.)
- 1x Burntec Dressing (4 in. x 4 in.)
- 5 x Burn Gel
- 1x MicroShield CPR Face Mask
- 5 x Antibiotic Ointment
- 3 x Sting & Bite Packet
- 3 x Salty Britches Chafing Ointment - Single Use
- 1x NAR® Emergency Survival Blanket
- 5 x Hand Sanitizer
- 2 pr Responder Nitrile Gloves, Large
- 1x Mini Duct Tape (2 in.)
- 1x Black Permanent Marker, Large
- 1x Bleeding Control Instruction Card
- 1x Triangular Bandage (40 x 40 x 56 in.)
- 1x Splinter Forceps, 4.5 in.
- 1x Trauma Shears, Large
- 1x Pelican Ruck R60 Case - OD Green
- 1x I-ROK Case Label



I-ROK XL

ITEM #	COLOR	OPTION
82-0049	Black	Basic
82-0051	OD Green	Basic
82-0053	Black	w/ BCD
82-0055	OD Green	w/ BCD
82-0070	Black	w/ CELOX™
82-0072	OD Green	w/ CELOX™



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277