

SOLO IFAK

Be prepared to save a life, maybe even your own

With NAR's SOLO IFAK (Individual First Aid Kit), you are equipped with combat proven treatments in an ultra-compact emergency solution. From the wilds of suburbia to the serenity of the mountains, never travel alone when you have the SOLO IFAK, this robust kit enables you to address your first aid needs. For injuries ranging from major bleeding trauma and splinting to minor cuts and wound care. Don't leave home without it.

IDEAL FOR

- Individual Sports
- Backpacking
- Hiking
- Hunting
- Camping
- Vehicles
- Adventure travel
- Surfing
- Motocross / ATV
- Mountain Biking
- Skiing
- Snowboarding



Don't leave home without it!

FEATURES/SPECIFICATIONS


- First aid for minor to severe injuries for single patient
- Ultra-compact & lightweight
- Optimized organization of contents
- Provide emergency treatment solutions until definitive care is available

Learn more about this and other **Products with a Mission®** at www.NARescue.com



SOLO IFAK

ITEM #	COLOR
82-0077	Black
82-0000	Coyote



DIMENSIONS

- H 7 in x W 4.5 in x D 2.75 in
- Weight: 154 lb



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

SOLO IFAK

KIT INCLUDES:



Learn more about this and other **Products with a Mission®** at www.NARescue.com



KIT CONTENTS

MAJOR TRAUMA

- 1x C-A-T® (Combat Application Tourniquet®) ORG
- 1x NAR® Responder ETD®, 4 in (Emergency Trauma Dressing)
- 1x Compressed Gauze, 4.5 in. x 4.1 yd
- 1x QuikClot® Bleeding Control Dressing
- 1x NAR® Emergency Survival Blanket - 56 in. x 96 in.

BANDAGES

- 1x BurnTec® Burn Dressing, 4 in x 4 in (10 cm x 10 cm)
- 10x Flexible Fabric Bandage, 0.75 in x 3 in
- 1x Wound Closure Strip (Pack of 10)
- 2x Gauze Pads, 4 in x 4 in (Pack of 2)
- 1x Clear Adhesive Dressing, 4 in x 5 in
- 1x Clear Adhesive Dressing, 2.4 in x 2.75 in
- 1x Silver Bandage, 2 in x 2 in
- 1x Silver Bandage, 4 in x 4.5 in
- 2x Petrolatum Gauze, 3 in x 18 in
- 1x NAR® Flat Tape

WOUND CARE

- 1x Irrigation Syringe, 10 cc

SPLINTING

- 1x Triangular Bandage, 40 in x 40 in x 56 in

EYE CARE

- 2x Polycarbonate Eye Shield (PES)

PERSONAL PROTECTION AND ACCESSORIES

- 1 pr Bear Claw® Nitrile Trauma Gloves-LG
- 1x NAR® Trauma Shears, Small
- 1x Permanent Marker, Blk, Small
- 1x Penlight
- 1x CPR Microshield Mask
- 1x Splinter Forceps, 4.5 in.
- 1x Luminous MED ID Patch



SOLO IFAK

ITEM #	COLOR
82-0077	Black
82-0000	Coyote



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277