

# BOA® IV CONSTRICTING BAND

## Increase the proficiency in gaining direct IV access...

When tasked with improving peripheral IV access, North American Rescue® engineers addressed both aspects of the challenge—how the product would work, and how easy it would be to work with. They met the challenge (and then some) with the BOA® Constricting Band, a now-critical component of successful IV cannulation.

How it works is revolutionary. Not since constricting bands were first conceived has the design been so fundamentally altered as it has been with the BOA. For improved central venous access every time, the BOA uses a rolling action to transfuse venous blood to the distal veins while still allowing arterial flow. (For maximum effect, the BOA is placed as high up on the extremity as allowable.) This technique distends the veins to make IV access quicker and more effective even in the most challenging patients.

How easy it is to work with is also revolutionary. The BOA's patented design features a Quick Click connector and two-finger Quick Release for added efficiency. No longer do healthcare professionals need to tie (and often re-tie) a constricting band. The BOA decreases "time on task" for treating combat casualties, civilian trauma patients, or medical patients with poor vascular access (for example, diabetic wound care).



### FEATURES/SPECIFICATIONS

- Quick, effective IV access even under challenging circumstances (including diabetic wound care)
- Revolutionary rolling action for transfusing venous blood to the distal veins while still allowing arterial flow
- Advanced design with Quick Click connector and two-finger Quick Release
- No tying required
- Rugged, vacuum-sealed packaging designed for sanitation and quick access (Single Use Only)
- US Patents: 5,304,202 and 5,607,448
- [Click here](#) to view the BOA Constricting IV Bands IFU

### DIMENSIONS

- Standard:*
- Device Length: L 7.5 in. to 18 in.
  - Packaged: L 4 in. x W 3 in. x D 0.9 in.
  - Weight: 1.2 oz
- XL:*
- Device Length: L 10 in. to 27 in.
  - Packaged: L 4 in. x W 3 in. x D 0.9 in.
  - Weight: 1.5 oz

Learn more about this and other **Products with a Mission®** at [www.NARescue.com](http://www.NARescue.com)



### BOA IV CONSTRICTING BAND

SIZE	ITEM #	SEL #	NSN #
Standard	30-0009	09MS-05-IVSA	6515-01-537-2611
XL	30-0071	09MS-05-IVSA	



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277