

# HEAVY DUTY TRAUMA XSHEARS

## XSHEAR AND XSHEAR MINI

The XSHEAR® is a heavy duty trauma shear designed to outperform traditional trauma shears in every way. The hardened 3mm thick stainless steel blades have a comfortable grip with slip resistant inner handles. The patented curved design for gentle approach when cutting near skin is hinged by a heavy duty center bolt that will not loosen over time. The razor sharp blades have a serrated lower blade to help cut through tough material with ease.

The XSHEAR® Mini Trauma Shear delivers all the benefits of the XSHEAR's premium cutting performance in a compact, lightweight design

Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



## FEATURES/SPECIFICATIONS

- Hardened stainless steel blades that are twice as thick as most trauma shears, sharpened to a razor edge
- BLACK TITANIUM coating on blades provides extra durability and sleek matte black appearance
- Patented blunt tip and curved design for gentler edge near the skin
- Soft touch inner handles for comfort and slip resistance
- Heavy duty center bolt that won't loosen over time
- **Mini Trauma Shears** deliver premium cutting performance in a compact, lightweight design

## DIMENSIONS

### Standard:

- Weight: 4 oz
- Length from longest point of handle to blade: 7.5 in.
- Width at widest point of handle: 3.5 in.
- Blade thickness approximately: 0.118 in.
- Cutting area length: 2 in.

### Mini:

- Weight: 1.55 oz
- Length from longest point of handle to blade: 6 in.
- Width at widest point of handle: 3 in.
- Blade thickness approximately: 0.115 in.
- Cutting area length: 2 in.

## XSHEAR® 7.5 IN.

ITEM #	COLOR
70-0290	Black/Black
70-0444	White/Rose Gold
70-0443	OD Green/Black
70-0291	Blue/Gray
70-0292	Orange/Black
70-0293	Pink/Black
70-0294	Red/Black
70-90415	Coyote/Black
70-424	Green/Black
70-0442	Purple/Black



## XSHEAR® MINI 6 IN.

ITEM #	COLOR
70-0435	Black/Black
70-0436	White/Rose Gold
70-0437	OD Green/Black
70-0438	Blue/Black
70-0441	Orange/Black
70-0439	Pink/Black
70-0440	Purple/Black



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277