

# TACTICAL CRICKIT®

## Specifically Engineered To Meet CoTCCC Requirements For A Surgical Airway Kit

The Tactical Crickit® is the result of an 18 month interface with Military Medical Professionals worldwide. Phase 1 of the project was to define the functional requirements of performing a surgical cricothyroidotomy in the tactical environment. Phase 2 was to initiate the development of a small, easy to use kit that focused on intuitive training behaviors. The result of this effort is the Tactical Crickit®.

Supplying the operator with all the necessary tools to conduct a proper cricothyroidotomy in a low cube and ruggedized package, the Tactical Crickit® is the ideal kit for performing a cricothyroidotomy in austere environments. For ease in both training and live fire situations, it's components are provided in the proper order of use. Once unrolled, the Tactical Crickit's 500-denier pouch provides a clean working surface for the operator, making it viable for use in low-light and battlefield environments. The Tactical Crickit also provides the user with a tracheal tube that is both the ideal size and length, along with a securing device to aid in secondary treatment and CASEVAC. To combat microbial contamination, the Tactical Crickit® is supplied with alcohol pads for immediate disinfection. The Tactical Crickit® was specifically engineered to meet CoTCCC requirements for a surgical airway kit.

**Learn more about this and  
other Products with a Mission®  
at [www.NARescue.com](http://www.NARescue.com)**



### FEATURES/SPECIFICATIONS

- The compact all-in-one cricothyroidotomy kit that contains all features required by CoTCCC for a surgical airway kit in the tactical setting
- Provides all of the components necessary to conduct a surgical cricothyroidotomy correctly
- The ideal length and size tracheal tube with securing device included
- Components configured in the sequence of use and ideal for low-light, battlefield environments.

### TACTICAL CRICKIT

ITEM #	NSN
10-0017	6515-01-540-7568

R<sub>X</sub> ONLY

### DIMENSIONS

- Kit Unrolled - H 6.5 in. x W 75 in. x D 1.25 in.
- Packaged - H 6.5 in. x W 2.5 in. x D 1.25 in.
- Weight: 3.7 oz



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

# TACTICAL CRICKIT®

**all** KITS INCLUDE:



**QTY 1**  
CricKit Nylon Pouch



**QTY 2**  
NAR Alcohol Prep Pads



**QTY 1**  
4" x 4" Gauze Pad 2 PK



**QTY 1**  
Protected #10 Scalpel



**QTY 1**  
NAR® Tracheal Hook



**QTY 1**  
Tube Securing Strap



**QTY 1**  
Cric. Tube, Cuffed, 6 mm ID



**QTY 1**  
Luer-Lock Syringe, 10cc

**KIT CONTENTS**

- 1 x Nylon Bag (Pouch- Tactical CricKit®)
  - 2 x NAR® Alcohol Prep Pads
  - 1 x Gauze pad (2 pack) - 4 in. x 4 in.
  - 1 x Protected #10 Scalpel
  - 1 x NAR® Tracheal Hook
- 1 x Tube Securing Strap
  - 1 x Cric. Tube with Obturator (Flexible, Cuffed, 6mm)
  - 1 x Luer-Lock Syringe - 10cc

Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



**TACTICAL CRICKIT**

ITEM #	NSN	Rx ONLY
10-0017	6515-01-540-7568	



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277