

NAR® TRIANGULAR BANDAGE

The NAR® Triangular Bandage is used for the bandaging of wounds, the immobilization of fractures and dislocations, and as a sling for the support of an injured part of the body.

Packaged for low cube space, the NAR® Triangular Bandage is made of muslin cloth and comes with two safety pins to aid in bandage configuration and securing. Easily fits into any first aid kit or pocket and is a staple for basic bandaging or splinting.

Learn more about this and other Products with a Mission® at www.NARescue.com



FEATURES/SPECIFICATIONS

- For the bandaging of wounds, the immobilization of fractures and dislocations, and as a sling for the support of an injured part of the body
- Packaged for low cube space
- Constructed of muslin cloth with two safety pins to aid in bandage configuration and securing
- Easily fits into any first aid kit or pocket

DIMENSIONS

- Packaged: 3.5 in. x 2 in. x .625 in.
- Opened: 37 in. x 37 in. x 52 in.
- Weight: 1.6 oz.

TRIANGULAR BANDAGE

ITEM #	NSN
30-0089	6510-00-201-1755



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277