

NAR® ELITE™ RESPONDER LITTER

Evacuation Litter for Incident Tactical Extraction

North American Rescue's Evacuation Litter for Incident Tactical Extraction (ELITE™) Responder Litter is a cost effective emergency carry litter designed for rapid deployment for the movement of casualties and is the ideal disposable litter for mass casualty incidents, disasters, active shooter events, as well as for everyday patient movements. Constructed of a water-resistant, woven polyethylene that is latex free and rated for up to 500 lb., the NAR ELITE™ Responder Litter features ten (10) carry handles with four (4) on each side and one (1) larger handle at either end of the litter to provide multiple carry points.



FEATURES/SPECIFICATIONS

- Cost effective emergency carry litter
- Features ten large carry handles with a load capacity tested to 500 lb
- Latex free, water-resistant

**Learn more about this and
other Products with a Mission®
at www.NARescue.com**



NAR® ELITE™ LITTER

ITEM #	NSN
60-0061	----

DIMENSIONS

- Litter Bed: L 76 in. x W 30 in.
- Packaged: L 11 in. x W 7.75 in. x D 2 in.
- Weight: 17 oz
- Payload: 500 lb



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277