

SQUAD BAG (GEN 2)

BAG ONLY

Learn more about this and other **Products with a Mission®** at www.NARescue.com



Designed for End-User Customization and Modernization

The Gen 2 Squad Bag is the latest innovation to the NAR® Squad Bag platform. The Gen 2 bag features about 3 in. of total width reduction and an added zipper pocket to expand or reduce internal cube space. The Gen 2 also features a 20% reduction in overall bag weight, with Laser Cut Laminate (LCL) Molle on both sides and front for additional carrying capabilities. Ultimately, end-user customization and modernization was at the forefront of the Gen 2 Squad Bag's design.

The Gen 2 can function as a fanny pack or over-the-shoulder bag, and is molle-attachable to other platforms. The bag is available in 500D Cordura Black or Coyote, and expandable up to a 56 in. waist.

The kit versions of the Squad Bag (Gen 2) feature off-the-shelf kit options to maximize your capabilities, including Basic Life Support (BLS), Advanced Life Support (ALS) and optional Hemostatic Dressing choices.



Multiple MOLLE attachment capabilities - additional hits sold separately.



Zipper Pocket to Expand or Reduce Internal Cube Space

FEATURES/SPECIFICATIONS

- Customizable, low-weight design with expandable cube-shaped interior
- Side straps can be adjusted to contour to body
- Durable 500D Nylon & Laser Cut Laminate Molle
- Can carry front fanny-style or backside
- Includes heavy-duty buckle hardware and self-repairing zippers with silent 550 cord pulls
- Sand & water resistant
- BAG ONLY. No components are included

SQUAD BAG (GEN 2)

ITEM #	COLOR	NSN
80-1793T	Black	----
80-1795T	Coyote	----

* Explore **Squad Kit (Gen 2)** options at www.NARescue.com.

DIMENSIONS

- BAG ONLY: H 6.25 in. x W 9.25 in. x D 6 in.
- Weight empty: 0.94 lb.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277