

IV POLE

ADJUSTABLE

Learn more about this and other **Products with a Mission®** at www.NARescue.com



Efficient and Reliable Solution...

North American Rescue's Adjustable IV Pole is exclusively designed to expand the treatment capabilities of our family of evacuation platforms by fitting the pre-molded holes of all NAR framed litter platforms including the Talon®, Raven®, T3, Medevac & OSL litters.

This essential, lightweight accessory is easy to transport and fully adjustable, from a minimum clearance of 18.5 in. expanding to a height of 32 in. when deployed.

Exclusively designed to fit the pre-molded holes of our framed litters and weighs less than 5 ounces.

NOTE: All accessories are sold separately. Litters are not included.

FEATURES/SPECIFICATIONS

- Essential, lightweight litter accessory designed to fit the pre-molded holes of all NAR framed litter platforms including the Talon®, Raven®, T3, Medevac & OSL litters
- Easy to transport and fully adjustable, from a minimum clearance of 18.5 in. expanding to a height of 32 in. when deployed

IV POLE

ITEM #	NSN
60-0008	6530-01-523-3314



DIMENSIONS

- Retracted: H 18.5 in. x W 8.86 in. x D 0.5 in.
- Fully Extended Height: 32 in.
- Weight: 4.9 oz.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277