

# PolarSkin™ Cooling Systems


INTRODUCING THE LATEST IN COOLING TECHNOLOGY FOR EFFECTIVE HUMAN THERMO-REGULATION

## THE LATEST COOLING TECHNOLOGY TO **BEAT the HEAT**

- Facilitates conductive, convective, and evaporative cooling
- Unique fiber construction allows for liquid absorptive capacity 350% greater than basis weight
- Protracted direct cooling of high blood flow regions
- Adheres to current National Athletic Trainers' Association and U. S. Ar t
- Soaked in ArcticAqua™ for faster and longer cooling without drying out



**PolarSkin Core includes a gel pack within the Arctic Mesh fabric that has been soaked in ArcticAqua to assist with cooling.**



**NORTH AMERICAN RESCUE®**  
www.NARescue.com • 888.689.6277

# PolarSkin Complete Care Backpack

*The Complete Care Backpack is a Comprehensive PolarSkin option designed to carry the Comprehensive Cooling System(s) into the field - with the ease and convenience of a Backpack style configuration.*

The PolarSkin Complete Care Backpack – Single Patient is a portable and highly insulated option that allows for the carry of One (1) PolarSkin Comprehensive Cooling System capable of treating a single patient in the field. The Cooler keeps ice from melting for 24 hours in 100°F heat. The PolarSkin™ Comprehensive Cooling System is a highly-effective portable cooling system for the treatment of heat-related illnesses. Consisting of two (2) PolarSkin™ Ice Sheets (Triple Pack) and one (1) PolarSkin™ C2E Ice Pack Bundle, the PolarSkin™ Comprehensive Cooling System is scientifically proven to promote rapid core temperature reduction.

**Single Patient Size Kit Contains:**

- 1 x PolarSkin Backpack Cooler, Black
- 1 x PolarSkin Comprehensive Cooling System containing:
  - 2 x PolarSkin Ice Sheets™ (Triple Pack)
  - 1 x PolarSkin C2E Ice Pack™ Bundle

**Double Patient Size Kit Contains:**

- 1 x PolarSkin Backpack Cooler, Black
- 2 x PolarSkin Comprehensive Cooling System each containing:
  - 2 x PolarSkin Ice Sheets™ (Triple Pack)
  - 1 x PolarSkin C2E Ice Pack™ Bundle

**DIMENSIONS**

*Single Packaged:*

- W 19.5 in. x L 15 in. x D 7.5 in.
- Weight: 12.8 lbs

*Double Packaged:*

- W 19.5 in. x L 15 in. x D 7.5 in.
- Weight: 22 lbs

**FEATURES/SPECIFICATIONS**

- Pack features Backpack Cooler and One (1) Comprehensive Cooling System for a portable, all-inclusive system capable of treating a single patient.
- Scientifically proven to reduce core temperature by .09°F/minute in 100°F heat when used according to protocol.
- The system maintains an internal temperature of ≤32°F for 17.5 hours in >85°F ambient temperatures when prepared to 0°F and stored in a PolarSkin Cooler with ice.
- Adheres to current National Athletic Trainers' Association and U.S. Army protocols for heat illness treatment

Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



**POLARSKIN COMPLETE CARE BACKPACK**

ITEM #	COLOR	SIZE	NSN
80-0543	Black	Single	----
80-0544	Black	Double	----



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277