

# USCG IFAK

## Designed To Solve Problems, Not Create Them

The trim USCG IFAK from North American Rescue® can be tucked into small spaces or carried about with ease until it is required for service. Built to the strict requirements of the United States Coast Guard, this improved first aid kit provides for safe first aid until a tactical healthcare professional arrives.

The USCG IFAK is a reliable solution. Made of the highest quality materials, our modular system can be worn by the operator as a drop-down leg rig, hooked directly to the duty belt, or attached to an assault vest or body armor with MOLLE-style connectors.

Learn more about this and  
other *Products with a Mission*®  
at [www.NARescue.com](http://www.NARescue.com)



## FEATURES/SPECIFICATIONS

- Durable 500D Nylon
- Sand & water resistant, self-repairing zipper with silent 550 cord pulls
- MOLLE-style Attachment System Compatible
- Removable Belt and Leg Straps
- Leg Strap with Quick-Grip™
- Sleek Design with no Front Pocket

## USCG IFAK

ITEM #	COLOR	NSN
80-0149	Black	6545-01-572-7248
80-0158	Coyote	6545-01-587-8513
85-0104	Resupply	----



## DIMENSIONS

- Case Closed: H 8.25 in. x W 71 in. x D 5.5 in.
- Weight: 2 lb 10 oz



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

# USCG IFAK

## KITS INCLUDE:



\*Resupply Kit, excludes Combat Casualty Card and Permanent Marker, and comes in a vacuum-sealed bag

## KIT CONTENTS

- 1x C-A-T® (Combat Application Tourniquet®) - Black
- 1x QuikClot® Combat Gauze®, Z-Fold Hemostatic Dressing - 3 in. x 4 yd.
- 1x NAR® S-Rolled Gauze - 4.5 in. x 4.1 yd.
- 1x ETD® Emergency Trauma Dressing - 6 in.
- 1x NAR® Pre-lubed Nasopharyngeal Airway - 28F (117 MM)
- 1x HyFin® Vent Chest Seal (Twin Pack)
- 1x Enhanced ARS® for Needle Decompression - 10 gauge x 3.25 in.
- 1x Polycarbonate (PES) Eye Shield
- 1x GeckoGrip® Rescue Tape (6 pack)
- 2 pr Bear Claw® Nitrile Trauma Gloves - Large
- 1x Combat Casualty Card\*\*
- 1x Permanent Marker - Large\*\*
- 1x Bag (USCG -IFAK)\*\*

\*\*Not included in Resupply Kit (85-0104)



## USCG IFAK

ITEM #	COLOR	NSN
80-0149	Black	6545-01-572-7248
80-0158	Coyote	6545-01-587-8513
85-0104	Resupply	----

Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277