

MEDIC ASSAULT RESCUE KIT

Middle hinge design allows contouring to the body

A Contoured Fit For Room Breaching And Clearing.

The Medic Assault Rescue Kit is our newest compact, low-profile trauma & rescue aid bag. This bag is unique because it is designed with a middle hinge, allowing a contoured fit for room breaching and clearing.

Two main bag compartments are accessible from each end of bag independently. The front of the bag has two large front pockets with full MOLLE on the exterior. Inside, there are four large mesh inner pouches that are secured to the bag interior with hook & loop attachment and provide organization via removable hook & loop labels for customization. The Basic Rescue kit comes with two Bleed pouches, one Splinting pouch, and one Airway pouch. The Advanced Rescue kit comes with two Bleed pouches, one IV pouch, and one Airway pouch. Additionally, the Rescue versions of the Medic Assault Kits come with carabiners, rope grabs, assorted ropes, and tubular webbing.

MOLLE on back and bottom of the exterior of the nylon bag provides additional carry capability. All silent zipper pulls. Removable backpack straps allow attachment directly to body armor if desired (requires additional straps/hardware – not included).

Available in six 500D Cordura colors, including MTC Alpine.

MEDIC ASSAULT RESCUE KIT



ITEM #	COLOR	OPTION	NSN
80-1737	Black	Basic	----
80-1738	Coyote	Basic	----
80-1740	OD Green	Basic	----
80-1741	Ranger Green	Basic	----
80-1739	Multicam	Basic	----
80-1736	MTC Alpine	Basic	----
80-1743*	Black	Advanced	----
80-1744*	Coyote	Advanced	----
80-1746*	OD Green	Advanced	----
80-1747*	Ranger Green	Advanced	----
80-1745*	Multicam	Advanced	----
80-1742*	MTC Alpine	Advanced	----

* Advanced kits are R_X ONLY

FEATURES/SPECIFICATIONS

- Slim Design and Low Empty Weight
- Middle hinge design to allow contouring to the body decreasing profile
- Large mesh internal bags
- Attachable to body armor
- Rescue versions come with carabiners, rope grabs, assorted ropes, and tubular webbing
- Available in six 500D Cordura colors, including MTC Alpine

DIMENSIONS

- L 22 in. x W 11.5 in. x D 5 in. (3 in. seam to seam)
- Basic Weight: 12.5 lbs.
- Advanced Weight: 13.8 lbs.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277



Learn more about this and other Products with a Mission® at www.NARescue.com



MEDIC ASSAULT RESCUE KIT

all KITS INCLUDE:



KIT CONTENTS

- 2 x C-A-T® Tourniquet - Black
 - 4 x NAR® Wound Packing Gauze - 3 in. x 5 yd.
 - 2 x NAR® Emergency Trauma Dressing (ETD®) - 6 in.
 - 2 x NAR® Flat Emergency Trauma Dressing (ETD®) - 4 in.
 - 1 x Abdominal/Stump Emergency Trauma Dressing (ETD®) - 12 in. x 12 in.
 - 2 x NAR® Nasopharyngeal Airway - 28FR - Pre-lubricated
 - 1 x Tactical Suction Device
 - 2 x HyFin® Vent Chest Seal Dressing (Twin Pack)
 - 2 x NAR® Emergency Survival Wrap
 - 1 x Ex-Tract Strap
 - 2 x SAM® Splint II
 - 2 x Triangular Bandage - OD Green
 - 2 x Polycarbonate Eye Shield (PES) - Orange
 - 1 x NAR® Trauma Shears, Black - Large - 7.25 in.
 - 1 x Trauma Shear Scissor Leash™ - Black
 - 2 x ROCKD Screw Lock Carabiner
 - 2 x ROCKO Auto-lock Carabiner
 - 2 x HMS Auto-lock Carabiner
 - 1 x Petzl® TiBloc Rope Grab / Ascender
 - 1 x F3 Anti-Panic Descent Device - Black
 - 1 x SPOC Rope Grab Pulley
 - 1 x Edelrid® HMPE Cord Blue Sling - 120 CM
 - 1 x Prusik VT Accessory Rope - 7 mm.
 - 1 x 50 ft. HAVOC Rope - 8 mm.
 - 1 x 75 ft. RIT 900 Rope - 6.8 mm.
 - 1 x Nylon Tubular Webbing, Black - 1 in.
 - 1 x Aviation Panel Kit - 35 in. x 29 in.
 - 2 x Combat Casualty Card
 - 2 x BLEED Nametape - Black/Gold
 - 1 x AIRWAY Nametape - Black/Gold
 - 1 x Permanent Marker, Black - Large
- plus, by kit type:**
- Basic**
- 1 x SPLINTING Nametape - Black/Gold
- Advanced**
- 2 x Quikclot® Z-Fold Combat Gauze - 3 in. x 4 yd.
 - 1 x Cyclone® Pocket Bag Valve Mask (BVM)
 - 1 x Bougie Aided Crick (BAC) Pack (Cricotrydotomy Kit)
 - 2 x Enhanced ARS™ Decompression Needle - 10 ga x 3.25 in.
 - 2 x Saline Lock Kit
 - 1 x BOA® IV Constricting Band
 - 1 x Intraosseous EZ-IO Manual Needle (Box of 6)
 - 1 x NAR® Mini Duct Tape - 2 in.
 - 1 x Halstead Straight Mosquito Forceps - 5 in.
 - 1 x Sharps Shuttle
 - 1 x Armadillo Medicine Case - Black
 - 1 x IV Nametape - Black/Gold



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277