

# EMERGENCY SURVIVAL WRAP

## The Ideal Solution For A Small, Lightweight Survival Wrap

The North American Rescue® Emergency Survival Wrap balances performance, durability and size. The wrap reflects radiant body heat and is constructed of a combination of mylar and a tear resistant polyethylene fabric material.

The wrap is wind and waterproof, further aiding in retaining warmth. The high-visibility, reversible orange and silver color helps rescuers find you during an emergency situation. Also available in reversible OD green & silver. The North American Rescue® Emergency Survival Wrap is an essential piece of equipment in a compact package to help maintain body heat.

**Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)**



## FEATURES/SPECIFICATIONS

- Reflects radiant body heat
- Lightweight, durable, laminated polyethylene
- Large 56 in. x 96 in. wrap for full body coverage
- Wind and waterproof
- Reversible OD green & silver color

## EMERGENCY SURVIVAL WRAP

ITEM #	COLOR
ZZ-0834	ORANGE
ZZ-1390	OD GREEN

## DIMENSIONS

- Packaged: L 5 in. x W 3.75 in.
- Orange Blanket: L 96 in. x W 60 in.
- OD Green Blanket: L 96 in. x W 56 in.
- Weight: 3.2 oz



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277