

# TRAUMA & FIRST AID BACKPACK KIT - CLASS B

## STATPACK® G3+ QUICKLOOK BAG FOR MOBILE RESPONSE

### NAR® TFAK Exceeds ANSI/ISEA Standards!

The **Trauma and First Aid Backpack Kit - Class B** is designed as a portable First Aid Kit with Class B trauma and first aid kit contents. It is available in three configurations: **Basic, Bleeding Control Dressing, or Celox™ Hemostatic Dressing**. Each kit **exceeds ANSI/ISEA Z308.1-2021** standards and includes additional bleeding control elements, such as the **C-A-T® Tourniquet, the fastest, safest, and most effective prehospital field tourniquet in the world**. Designed to treat the most common workplace injuries while also controlling major bleeding, the Class B kit is stocked with a large quantity of supplies for larger or higher-risk workplaces.

The **Statpack G3+ Quicklook bag** allows versatility and the ability to be a **complete mobile solution for workplace and community environments**. When paired with a **ZOLL® Plus AED** and **NAR's BurnTec Emergency Burn Kit** (each sold separately), you have the ability to provide emergency care from a minor injury to cardiac arrest.

Be prepared and protected when rendering aid wherever you are. **Biological hazards** like bodily fluids may sometimes require more protection than gloves that are typically found in first aid and bleeding control kits. **NAR's Bloodborne Pathogen Personal Protection Equipment (PPE) Kit** is a great addition to any kit for extra protection to have available when needed.

*For workplace requirements OSHA's Bloodborne Pathogens standard can be found in Title 29 of the Code of Federal Regulations at 29 CFR 1910.1030. Part of this standard requires employers to provide personal protective equipment (PPE), such as gloves, gowns, eye protection, and masks. All of which can be found in NAR's Bloodborne Pathogen Kit.*

### FEATURES/SPECIFICATIONS

- **Class B** for larger workplaces or higher-risk environments
- Configurations: **NAR® Compressed Gauze, QuikClot® Bleeding Control Dressing, or Celox™ Hemostatic Dressing**
- **Exceeds ANSI/ISEA Standards**
- Able to **treat the most common workplace injuries** with the added benefit of being able to **control major bleeding** with the C-A-T® Tourniquet, the fastest, safest, most effective prehospital field tourniquet in the world

- **Statpack® G3+ Quicklook Bag** for mobile response
- Compatible with **ZOLL® Plus AED & BurnTec Kit** (sold separately)
- **Bloodborne Pathogen Kit** added for extra protection

*Zoll® AED  
compatible.  
Sold separately*



**Learn more about this and  
other Products with a Mission®  
at [www.NARescue.com](http://www.NARescue.com)**



### TRAUMA & FIRST AID BACKPACK KIT - CLASS B

ITEM #	OPTION
80-1857	NAR® Compressed Gauze
80-1861	QuikClot® Bleeding Control Dressing
80-1865	Celox™ First Aid Gauze

### DIMENSIONS

- H 19 in. x W 17 in. x D 9 in.
- Weight: 10.5 lbs




**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

# TRAUMA & FIRST AID BACKPACK KIT - CLASS B



**PLUS+**  
one of these:



**Removable  
internal  
pouches!**



## TRAUMA & FIRST AID BACKPACK KIT - CLASS B

ITEM #	OPTION
80-1857	NAR® Compressed Gauze
80-1861	QuikClot® Bleeding Control Dressing
80-1865	Celox™ First Aid Gauze

### KIT CONTENTS

- 1x C-A-T® Tourniquet – Orange
- 1x NAR® Responder Compressed Gauze
- 1x NAR® Flat Responder ETD® - 4 in.
- 1x CPR Microshield Mask
- 1x HyFin® Vent Compact Chest Seal Twin Pack
- 1x Emergency Survival Blanket
- 2x Triangular Bandage
- 1x SAM Universal Malleable Splint
- 1x NAR® Elastic Wrap Bandage – 4 in.
- 2x NAR® Elastic Wrap Bandage – 2 in.
- 4x Gauze Pad (2 pack) – 4 in. x 4 in.
- 4x Gauze pad – 5 in. x 9 in.
- 2x BurnTec® Dressing – 4 in. x 4 in.
- 25x Safetec Burn Gel
- 50x Antiseptic Towelette
- 25x Antibiotic Ointment
- 20x Hand Sanitizer Packet
- 2x Instant Cold Pack – 6 in. x 9 in.
- 1x Eye Wash Solution – 1 oz.
- 2x Sterile Oval Eye Pad
- 2x Surgical Tape – 1 in.
- 50x Fabric Bandage – 1 in. x 3 in.
- 4pr Responder Gloves – Large
- 1x NAR® Responder Shears – Large
- 1x Bloodborne Pathogen PPE Kit:
  - 1pr Responder Gloves – Large
  - 1pr Responder Gloves – Extra Large
  - 2x Antiseptic Towelettes
  - 1x Biohazard Bag
  - 1x Impervious Gown – Universal Size
  - 1x Surgical Mask with Visor
  - 2x Sani-Hands® Sanitizing Wipes
  - 1x NAR® Sharps Container
- 1x Black Permanent Marker – Large
- 1x First Aid Pocket Reference Guide
- 1x Bleeding Control Instructions
- 1x Class A Kit Contents Inventory Card
- 1x Inspection Card

### Plus one of these:

- 1x NAR® Responder Compressed Gauze
- 1x QuikClot® Bleeding Control Dressing
- 1x Celox™ First Aid Hemostatic Dressing

Learn more about this and  
other Products with a Mission®  
at [www.NARescue.com](http://www.NARescue.com)



**TFAK Backpacks exceed ANSI/ISEA Z308.1-2021 standards and  
includes additional bleeding control elements!**



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277