

KOBUS ROLLUP LITTER STRAPS

TOW STRAP | SHOULDER STRAP (2PC) | FOOT STRAP (2PC)

Kobus Rollup Litter Tow Strap

Kobus Tow Strap with 2 handles & 2 carabiners with adjustable length. Tow Strap with quick release carabiners, adjustable for different pulling setups: single person, two people side by side, or one forward/one aft (requires two straps for ski configuration). Enables efficient use of torso, waist, and arms to reduce fatigue in both single and dual user modes. Includes a MOLLE/PALS compatible pouch.



Kobus Rollup Litter Shoulder Straps (2pc)

Kobus Shoulder Straps for patient positioning and fixing with adjustable length. Designed to secure patients at the shoulders and prevent sliding toward the head of the RollUp Litter, especially during high angle rescues. These adjustable straps replace time-consuming rope alternatives, can be quickly attached as needed. Includes MOLLE/PALS compatible pouch.



Kobus Rollup Litter Foot Strap (2pc)

Kobus Foot Strap Set (FSS) for patient positioning and fixing with adjustable length. Eliminates having to create straps out of rope/line which is time consuming and not easily adjusted. The FSS can be secured in less than 3 minutes and when affixed, easily positions and adjusts to the patient minimizing longitudinal movement towards the foot end of the RollUp Litter, an absolute necessity for shorter individuals. The FSS includes a mesh pouch which is easily stored with other accessories in/on the RollUp carry case. MOLLE/PALS compatible.



FEATURES/SPECIFICATIONS

- Kobus RollUp compatible (sold separately)
- Adjustable Length
- **Tow Strap**
 - » Includes 2 Handles & 2 Carabiners
 - » Quick Release Carabiners
 - » Includes MOLLE/PALS Compatible Pouch
- **Shoulder straps**
 - » For patient positioning and fixing
 - » Quick attachment
 - » Includes MOLLE/PALS Compatible Pouch
- **Foot Strap Set (FSS)**
 - » For Patient Position and Fixing
 - » Quick Installation
 - » MOLLE/PALS Compatible

DIMENSIONS

Tow Strap:

- Packed Dimensions: H 3 in. x L 8 in. x W 6 in.
- Weight: 1 lb.

Shoulder Strap:

- Packed Dimensions: H 4 in. x L 9 in. x W 7 in.
- Weight: 1 lb.

Foot Strap:

- Packed Dimensions: H 3 in. x L 8 in. x W 6 in.
- Weight: 1 lb.

KOBUS STRAPS

ITEM #	COLOR	ITEM	NSN
60-0114	Black	Tow Strap	----
60-0112	Green/Black	Shoulder Straps	----
60-0113	Red/Yellow	Shoulder Straps	----
60-0100	Green/Black	Foot Straps	----
60-0101	Red/Yellow	Foot Straps	----

* More Kobus products available at www.NARescue.com

Learn more about this and other Products with a Mission® at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277