

KOBUS ENCLOSED ROLLUP RESCUE STRETCHER

The optimal solution for rescue in confined spaces, heights, depths and horizontal rescues.

The patented RollUP system combines the advantages of a roll-up stretcher with the rigidity of a basket carrier.

The Kobus Enclosed RollUP Rescue Stretcher is designed after the Kobus RollUP 4000 series but with a slimmer form. It can be used for hoisting by helicopter, hoisting in caves or industrial confined spaces.

Once the patient is secured, the stretcher becomes rigid. Integrated struts prevent the stretcher from buckling and prevents pressure being exerted on the patient during the rescue.

The integrated struts can be removed before use if necessary - especially in tight spaces. All components are designed for chemical resistance and easy cleaning (decontamination or disinfection).

For convenience, the stretcher can be rolled up and stored in an included backpack.

FEATURES/SPECIFICATIONS

- Based on Kobus 4000 Series with slimmer form
- Hoistable. Easy transition between vertical and horizontal hoist
- Fits in confined spaces
- Integrated struts prevent buckling or pressure on the patient
- Material is chemical resistant
- Easily cleaned
- Roll up and store in an included backpack.

KOBUS ROLLUP RESCUE STRETCHER

ITEM #	COLOR	WIDTH	NSN
60-0105	Green	27.5 in.	----
60-0107	Green	36.2 in.	----
60-0109	Red	27.5 in.	----
60-0111	Red	36.2 in.	----

* More Kobus products available at www.NARescue.com

Learn more about this and other Products with a Mission® at www.NARescue.com



DIMENSIONS

27.5 in"

- Product Dimensions: L 96 in. x W 27.5 in.
- Weight: 11 lb..


36.2 in"

- Product Dimensions: L 97.6 in. x W 36.2 in.
- Weight: 18.3 lb.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

KOBUS ENCLOSED ROLLUP RESCUE STRETCHER



KOBUS ROLLUP RESCUE STRETCHER

ITEM #	COLOR	WIDTH	NSN
60-0105	OD Green	27.5 in.	----
60-0107	OD Green	36.2 in.	----
60-0109	Red	27.5 in.	----
60-0111	Red	36.2 in.	----

* More Kobus products available at www.NARescue.com

**Learn more about this and
other Products with a Mission®
at www.NARescue.com**



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277