

# KOBUS ROLLUP RESCUE STRETCHER

Learn more about this and  
other **Products with a Mission®**  
at [www.NARescue.com](http://www.NARescue.com)



**The optimal solution for rescue in confined spaces,  
heights, depths and horizontal rescues.**

*The patented RollUP system combines the advantages of a roll-up stretcher with the  
rigidness of a basket carrier.*

The Kobus RollUP Rescue Stretcher is designed after the Kobus RollUP  
4000 series but with a slimmer form. It can be used for hoisting by  
helicopter, hoisting in caves or industrial confined spaces.

Once the patient is secured, the stretcher becomes rigid. Integrated struts prevent the  
stretcher from buckling and prevents pressure being exerted on the patient during the rescue.

The integrated struts can be removed before use if necessary -  
especially in tight spaces. All components are designed for chemical  
resistance and easy cleaning (decontamination or disinfection).

For convenience, the stretcher can be rolled up and stored in an included backpack.

## FEATURES/SPECIFICATIONS

- Based on Kobus 4000 Series with slimmer form
- Hoistable. Easy transition between vertical and horizontal hoist
- Fits in confined spaces
- Integrated struts prevent buckling or pressure on the patient
- Material is chemical resistant
- Easily cleaned
- Roll up and store in an included backpack

## KOBUS ROLLUP RESCUE STRETCHER

ITEM #	COLOR	WIDTH	NSN
60-0104	OD Green	27.5 in.	----
60-0106	OD Green	36.2 in.	----
60-0108	Red	27.5 in.	----
60-0110	Red	36.2 in.	----

\* More Kobus products available at [www.NARescue.com](http://www.NARescue.com)



**Integrated struts on all  
Kobus Rescue Stretchers  
prevent buckling or  
pressure on the patient**



## DIMENSIONS

27.5 in"

- Package Dimensions: H: 10.6 in. x W: 28 in.
- Product Dimensions: L 96 in. x W 27.5 in.
- Weight: 121 lb.


36.2 in"

- Package Dimensions: H: 10.6 in. x W: 36.2 in.
- Product Dimensions: L 97.6 in. x W 36.2 in.
- Weight: 174 lb.



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

# KOBUS ROLLUP RESCUE STRETCHER



## KOBUS ROLLUP RESCUE STRETCHER

ITEM #	COLOR	WIDTH	NSN
60-0104	Green	27.5 in.	----
60-0106	Green	36.2 in.	----
60-0108	Red	27.5 in.	----
60-0110	Red	36.2 in.	----

\* More Kobus products available at [www.NARescue.com](http://www.NARescue.com)

**Learn more about this and  
other Products with a Mission®  
at [www.NARescue.com](http://www.NARescue.com)**



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277