

KOBUS ROLLUP LITTER INFLATABLE PILLOW

Kobus presents a compact and durable inflatable pillow which provides stabilization and positioning assistance of the patients head when properly applied. The pillow is added during the RollUp assembly process with no additional time required. It is then inflated when needed, otherwise, it stays out of the way. Inflation time (oral life jacket style inflation) is <15 seconds. As it is inflated the pillow provides cushion to the back of the head while grasping the sides of the head for stabilization. A forehead strap may be applied as needed. Includes a MOLLE/PALS compatible pouch.

Learn more about this and other Products with a Mission® at www.NARescue.com



FEATURES/SPECIFICATIONS

- Inflatable Pillow
- Kobus RollUp compatible (sold separately)
- For Patient Stabilization and Positioning
- Inflate Only When Needed
- Quick Inflation
- Includes MOLLE/PALS Compatible Pouch

DIMENSIONS

- Packed Dimensions: H 4 in. x L 9 in. x W 7 in.
- Weight: 1 lb.

KOBUS INFLATABLE PILLOW

ITEM #	COLOR	NSN
60-0102	Green	----
60-0103	Red	----

* More Kobus products available at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277