

KOBUS ROLLUP LITTER INFLATABLE PILLOW

Kobus presents a compact and durable inflatable pillow which provides stabilization and positioning assistance of the patients head when properly applied. The pillow is added during the RollUp assembly process with no additional time required. It is then inflated when needed, otherwise, it stays out of the way. Inflation time (oral life jacket style inflation) is <15 seconds. As it is inflated the pillow provides cushion to the back of the head while grasping the sides of the head for stabilization. A forehead strap may be applied as needed.

Kohlbrat and Bunz Corporation (Kobus) has a 10 year Manufacturer's Defect Warranty for Kobus products. For more information contact NAR® Customer Solutions.

Learn more about this and other Products with a Mission® at www.NARescue.com



FEATURES/SPECIFICATIONS

- Inflatable Pillow
- Kobus RollUp compatible (sold separately)
- For Patient Stabilization and Positioning
- Inflate Only When Needed
- Quick Inflation

DIMENSIONS

- Packed Dimensions: H 4 in. x L 9 in. x W 7 in.
- Weight: 1 lb.

KOBUS INFLATABLE PILLOW

| ITEM # | COLOR | NSN |
|---------|-------|------|
| 60-0102 | Green | ---- |
| 60-0103 | Red | ---- |

* More Kobus products available at www.NARescue.com

Compatible with Kobus Rollup Litters (sold separately)

