

# DERMAQUIK®


## SPORT & EXTREME (3PK)

Learn more about this and  
other **Products with a Mission®**  
at [www.NARescue.com](http://www.NARescue.com)

The DermaQuik® Sport & Extreme 3 Pack is a perfect addition to First Aid Kits for emergency laceration closure for the wilderness, work-site, sports field, campground, or wherever first aid kits are needed or access to facilities is limited. The Sport & Extreme device utilizes a more aggressive adhesive intended to adhere in harsh or austere conditions.

This adhesive is a breathable, non-sensitizing, long-duration adhesive intended for wear up to 10 days. The adhesive is acrylic based and pressure sensitive, so be sure to firmly massage the adhesive into the skin at application. This adhesive is not suitable for use on fragile or geriatric skin.

Each multi-pack includes 3 DermaQuik® Sport & Extreme devices, each closing 1 in. of wound (or use multiples for long wounds).



**Close Wounds  
in Seconds!**



### FEATURES/SPECIFICATIONS

- A non-invasive, needle-free alternative to sutures and glues for small cuts and lacerations
- Close wounds in seconds
- Potentially avoid Emergency Department (ED) visits, wait times and expensive hospital bills
- Close skin with one hand, in case no one else is around
- No mess of glue runoff, no trauma of suture puncture, no micro-needles
- Closes THREE 1 in. cuts or ONE 3 in. cut

### DERMAQUIK®

ITEM #

ZZ-1975

### DIMENSIONS

- L 6.75 in. x W 4.25 in. x D 0.2 in.
- Weight: 1 oz.



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277