

Semi-Recessed Metal Wall Case

An Ideal Solution To Safely Store & Protect

North American Rescue's Semi-Recessed Metal Wall Case provides low profile, high visibility access to critical medical supplies for use by initial first responders and bystanders in public areas. It is an ideal solution to safely store & protect Automated External Defibrillators (AEDs), Class A & B Trauma and First Aid Kits, and Public Access Bleeding Control Kits.

The cold steel cabinet with a white finish is ADA compliant and features concealed hinges, magnetic closure & recessed handle designed to open with one quick pull, without the need to twist and pull. The wall-mounted design fastens permanently inside the wall with a low profile 3.25 in. trim. Features include highly visible lettering (AED, FIRST AID KIT, or BLEEDING CONTROL KIT) for fast identification, and a clear plexiglass window for viewing cabinet contents. The Bleeding Control Kit cabinet also includes convenient Instructions for Use posted on the interior back panel. The optional 80db keyed alarm system, housed in a durable internal closure, offers security for stored medical assets, alerting all around deployment or tampering. Cases are water resistant. To best protect your critical life saving products, cases are best mounted under cover and out of the direct elements.

**Learn more about this and
other Products with a Mission®
at www.NARescue.com**



FEATURES/SPECIFICATIONS

- White powder coated cold rolled steel cabinet with a recessed handle, concealed hinges, and magnetic closure designed to safely store medical kits and equipment in public areas
- Provides both physical and environmental protection for medical equipment
- Easily mounts to any wall
- ADA Compliant
- Built for easy access - one quick pull to open, no twist and pull
- Designed to house
 - » All AED models from major manufacturers,
 - » NAR Class A & B Trauma & First Aid Kits
 - » NAR Public Access Bleeding Control Kits from North American Rescue
- Features highly visible lettering for fast identification:
 - » AED
 - » FIRST AID KIT
 - » BLEEDING CONTROL KIT
- Optional 80db keyed alarm system

SEMI-RECESSED METAL WALL CASE

OPTION	ITEM #
AED	80-1775
AED w/ Alarm	80-1776
First Aid Kit	80-1777
First Aid Kit w/ Alarm	80-1778
Bleeding Control Kit	80-1779
Bleeding Control Kit w/ Alarm	80-1780

* Explore Public Access Bleeding Control Stations and other cabinet options at www.NARescue.com

DIMENSIONS

- External Lip: H 16 in. x L 17 in. x D 3.25 in.
- Recessed Exterior: H 14 in. x W 15 in. x D 7.25 in.
- Rough wall opening: H 15 in. x L 16 in. x D 4.25 in.
- Weight: 13 lbs.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277