

FOX EYE SHIELD *with included Garter*



Convex Rigid Eye Shield...

The Rigid Eye Shield meets current military and civilian pre-hospital guidelines for the treatment of traumatic eye injury in the field.

This convex rigid eye shield allows the user to protect the injured eye and blocks external pressure from being applied directly to the injury, preventing further damage. It is made of aluminum for rigidity, yet contains many apertures to allow for air circulation and to reduce weight.

The garter (included) increases patient comfort as it pads the edges of the device upon application.

Features:

- Combination Eye Garter Shield Cover and Fox Shield; includes E5692 Aluminum Fox Eye Shield and One E569W Eye Garter Shield Cover White Only
- Lightweight aluminum
- Perforated to allow air circulation
- Accommodates either eye
- Meets CoTCCC treatment protocols for penetrating eye trauma
- Garter included for patient comfort

Dimensions:

- Dimensions:
2.938 in. x 2.312 in.
- Individual Weight:
0.2 oz (5 grams)

Fox Eye Shield

Option	Qty	Item #	NSN#
w/ Garter	Single	30-0055	
w/ Garter	Box of 50	30-5550	6510-01-598-1877

CE

Allows the user to protect the injured eye and blocks external pressure.

PRODUCTS WITH A
MISSION[®]



NORTH AMERICAN RESCUE[®]

www.NARescue.com • 888.689.6277