

LEASHED TRAUMA KIT

Lifesaving Products Within Arm's Reach

Introducing NAR's Leashed Trauma Kit featuring the Combat Application Tourniquet® (C-A-T®). Designed for those in environments where keeping lifesaving products within arms reach is paramount.

Our all new leashed system allows attachment with MOLLE backing to multiple platforms. The high viz orange quick pull tab gives the user easy one-handed access to an interior shock-cord panel made for life saving equipment. The panel is attached to the bag with a bungee type cord to keep lifesaving components within your workspace.

**Learn more about this and
other Products with a Mission®
at www.NARescue.com**



LEASHED TRAUMA KIT

ITEM #

80-1933

DIMENSIONS

- Packaged: L 7 in. x W 3 in. x D 2.5 in.
- Weight: 9 oz.

FEATURES/SPECIFICATIONS

- Protective hood with quick latch release helps protect contents and includes sewn in embroidered medical cross identifier
- Quick pull tab for easy one-handed deployment of shock-cord panel
- Allows attachment via MOLLE backing
- Shock-cord inner panel adjusts to secure internal components & is easily removable for rapid access



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

LEASHED TRAUMA KIT

all KITS INCLUDE:



LEASHED TRAUMA KIT

ITEM #

80-1933

Learn more about this and other Products with a Mission® at www.NARescue.com



KIT CONTENTS

- 1x Combat Application Tourniquet® (C-A-T®) - Orange
- 1x NAR® Mini Responder Emergency Trauma Dressing (ETD®) - 4 in.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277