

Stop the Lethal Diamond of Death with Warm Blood!

# INSTANT. INTUITIVE. PORTABLE.

## °MEQU °M WARMER SYSTEM

Portable Blood & IV-Fluid Warmer

**EASY SETUP:** User-friendly design;  
Standard Luer Locks

**LIGHTWEIGHT:** Battery with warmer only 1.8 lbs

**INSTANT WARMING:** Warms fluids in less than  
10 seconds

**MEETS COTCCC GUIDELINES:** Warms fluids to  
body temperature at flow rates up to 150 ml/min

**LARGE WARMING CAPACITY:** Warms up to 2L of  
cold fluids on one fully charged battery

**LOW MAINTENANCE:** Wipe clean and recharge  
after use

**FIELD-PROVEN:** In use since 2018 in both  
Military, HEMS, EMS and Hospitals

**DUST & SPLASH PROOF:** IP54 battery and  
warmer

**ENTIRE CONTINUUM OF CARE:** °M Warmer can  
stay with the patient for up to 72 hours

**ADHESIVE PATCH:** For secure attachment close  
to the infusion site

**INTEROPERABLE:** Seamlessly supporting  
continuity of care across the entire patient journey



### °M WARMER

Sterilized single patient use device that  
heats fluid/blood

**Dimensions:** 5 cm x 10 cm x 2 cm  
**Weight:** 0.24 lbs



### POWER PACK PLUS

Rechargeable 21.6V battery that supplies  
power to the °M Warmer

**Dimensions:** 9 cm x 18 cm x 3.5 cm  
**Weight:** 1.56 lbs



REV100925

### OPERATING PARAMETERS

<b>Output temp. range</b>	39°C/102.2°F (+/- 3°C/5.4°F at flowrates from 2-150ml/min)
<b>Max fluid warming volume on one fully charged Power Pack Plus (output temperature of min. 36°C/96.8°F)</b>	Input temp. 5°C/41°F: 2L Input temp. 21°C/69.8°F: 4L
<b>Priming volume</b>	3.5 ml (including tubing and Luer Locks)
<b>Maximum inlet pressure</b>	300 mmHg

DESCRIPTION	ITEM #	NSN
Starter Kit Plus	35-0048	6515-01-723-1196
Power Pack Plus	35-0049	6515-01-723-1180
°M Warmer - 5 pc	35-0050	6515-01-723-0591
DC Charger	35-0051	--
AC Charger	35-0052	6515-01-723-1257
°M Warmer System Pouch	35-0053	--



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

Discover the °M Blood &  
IV-Fluid Warmer & more  
Products with a Mission®  
at [www.NARescue.com](http://www.NARescue.com)

