

TRAUMA & FIRST AID KITS (TFAK) HARD CASE - CLASS A

The workplace Hard Case First Aid Kits that exceed ANSI/ISEA Z308.1-2021 for Class A

The TFAK (Trauma and First Aid Kit) - Class A from North American Rescue® is designed as a portable First Aid Kit for the workplace that exceeds ANSI/ISEA Z308.1-2021 Standard for Class A kits. Our kits exceed the ANSI and OSHA standards by adding additional bleeding control elements.

Packaged in a rigid carrying case that is easily transported to the patient's side, this kit offers the ability to treat the most common workplace injuries with the added benefit of being able to control major bleeding with the safest and most effective prehospital tourniquet on the market today.

FEATURES/SPECIFICATIONS

- Workplace first aid kits that exceed ANSI/ISEA Z308.1-2021 Standard for Class A kits
- Able to treat the most common workplace injuries with the added benefit of being able to control major bleeding with the safest and most effective Tourniquet on the market
- Designed so the kit can be easily carried to the patient's side

TRAUMA & FIRST AID KITS HARD CASE - CLASS A

ITEM #	COLOR	OPTION
80-1033	RED	w/Compressed Gauze
80-1046	RED	w/Bleeding Control Dressing
80-1481	RED	w/CELOX™ gauze
80-1914	ORG	w/Compressed Gauze
80-1915	ORG	w/Bleeding Control Dressing
80-1916	ORG	w/CELOX™ gauze
80-1917	YLW	w/Compressed Gauze
80-1918	YLW	w/Bleeding Control Dressing
80-1919	YLW	w/CELOX™ gauze

Modular system featuring a rugged storage organizer in the lid and removable nylon pouches with easy-to-read treatment labels...

DIMENSIONS

- H 6 in. x W 13.5 in. x D 11.5 in.
- Weight: 6 lbs.

Learn more about this and other **Products with a Mission®** at www.NARescue.com



ANSI

Type IV ANSI Containers

... With the added benefit of major bleeding control with the safest, most effective tourniquet on the market today!





NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

TRAUMA & FIRST AID KITS (TFAK) HARD CASE - CLASS A



**PLUS*
one of these:**



Learn more about this and other Products with a Mission® at www.NARescue.com



**Exceeds ANSI/ISEA Z308.1-2021
Standard for Class A Kits**



KIT CONTENTS

- 1 x C-A-T® Tourniquet, Orange
 - 1 x HyFin® Vent Chest Seal Twin Pack
 - 1 x NAR® Responder Compressed Gauze
 - 1 x 4 in. NAR® Responder ETD®
 - 1 x 4 in. Elastic Wrap Bandage
 - 1 x 2 in. Elastic Wrap Bandage
 - 1 x Triangular Bandage
 - 1 x SAM® Splint
 - 2 x Gauze, Pad - 4 in. x 4 in. (Pack of 2)
 - 2 x Gauze, Pad - 5 in. x 9 in.
 - 16 x Flexible Fabric Bandages, 1 in. x 3 in.
 - 1 x BurnTec® Dressing - 4 in. x 4 in.
 - 1 x CPR Microshield Mask
 - 10 x Antiseptic Towelette
 - 10 x Antiseptic Ointment
 - 1 x Roll 1 in. Surgical Tape
 - 1 x Instant Cold Pack - 6 in. x 9 in.
 - 1 x Eye Wash Solution, 1 oz
 - 2 x Sterile Oval Eye Pad
 - 2 pr Responder Nitrile Gloves - Large
 - 1 pr Responder Trauma Shears - Large
 - 10 x Hand Sanitizer
 - 10 x Burn Gel
 - 1 x NAR® Survival Blanket
 - 1 x Black Permanent Marker, Large
 - 1 x First Aid Pocket Guide
 - 1 x Just In Time Bleeding Control Instructions
 - 1 x Trauma & First Aid Kit Contents Inventory Card
 - 1 x Inspection Card
- Plus one of these:**
- 1 x NAR® Responder Compressed Gauze
 - 1 x QuikClot® Bleeding Control Dressing
 - 1 x Celox™ First Aid Hemostatic Dressing



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277