

RANGE TRAUMA KIT (HARD CASE)

Hardcase emergency first aid kit for treating injuries on the range

Treatment capabilities include penetrating & blast injuries, chest wounds, burns, eye injuries, minor lacerations, sprains and fractures and many other common injuries seen on training and shooting ranges.

Bright orange, water-tight case maintains the integrity of your medical products when you need them most.

Kit contents are arranged in a modular system featuring a rugged storage organizer in the lid and removable nylon pouches with easy-to-read treatment labels in the main compartment.



FEATURES/SPECIFICATIONS

- Brightly colored, water-tight case
- Modular system featuring a rugged storage organizer in the lid and removable nylon pouches with easy-to-read treatment labels



Learn more about this and other Products with a Mission® at www.NARescue.com



RANGE TRAUMA KIT (HARD CASE)

ITEM #	COLOR	OPTION
80-1901	Orange	Basic
80-1902	Red	Basic
80-1903	Yellow	Basic
80-1907	Orange	w/ Bleeding Control Dressing
80-1908	Red	w/ Bleeding Control Dressing
80-1909	Yellow	w/ Bleeding Control Dressing
80-1904	Orange	w/ CELOX™
80-1905	Red	w/ CELOX™
80-1906	Yellow	w/ CELOX™

DIMENSIONS

- L 9.5 in. x W 13.5 in. x H 6 in.
- Weight: 8 lb 12 oz




NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

RANGE TRAUMA KIT (HARD CASE)

all KITS INCLUDE:



+plus, by kit type:



Learn more about this and other Products with a Mission® at www.NARescue.com



KIT CONTENTS

- 2 x C-A-T® Tourniquet (Orange)
- 2 x NAR® S-Rolled Gauze
- 2 x 6 in. Emergency Trauma Dressing (ETD™)
- 1 x HyFin® Vent Chest Seal Twin Pack
- 5 x 2-Pack Sterile Gauze Pads, 4 in. x 4 in.
- 2 x Polycarbonate Eye Shield
- 2 x Elastic Wrap, 4 in.
- 2 x Triangular Bandage
- 2 x SAM® Splint II
- 1 x Surgical Tape, 2 in.
- 50 x Flexible Adhesive Strip Bandage, 0.75 in. x 3 in.
- 1 x BurnTec® Burn Dressing, 5 in. x 5 in.
- 1 x Eye Wash, 4 oz
- 1 x Forceps, Splinter, 4.5 in.
- 1 x Responder Trauma Shears, 7.25 in.
- 4 x pair Responder Nitrile Gloves, Lg.
- 1 x CPR Mask
- 1 x Marker, Permanent

+plus, by kit type:

- 2 x NAR® S-Rolled Gauze
- or -
- 2 x Combat Gauze® Hemostatic Dressing
- or -
- 2 x Celox™ First Aid Z-Fold Hemostatic Dressing



RANGE TRAUMA KIT (HARD CASE)

ITEM #	COLOR	OPTION
80-1901	Orange	Basic
80-1902	Red	Basic
80-1903	Yellow	Basic
80-1907	Orange	w/ Bleeding Control Dressing
80-1908	Red	w/ Bleeding Control Dressing
80-1909	Yellow	w/ Bleeding Control Dressing
80-1904	Orange	w/ CELOX™
80-1905	Red	w/ CELOX™
80-1906	Yellow	w/ CELOX™



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277