

TRAUMA & FIRST AID KITS (TFAK) - CLASS A

The workplace First Aid Kits that exceed ANSI/ISEA Z308.1-2021 for Class A...

The TFAK (Trauma and First Aid Kit) - Class A from North American Rescue® is designed as a portable First Aid Kit for the workplace that exceeds ANSI/ISEA Z308.1-2021 Standard for Class A kits. Our kits exceed the ANSI and OSHA standards by adding additional bleeding control elements. Packaged in a convenient zippered nylon bag that is easily transported to the patient's side, this kit offers the ability to treat the most common workplace injuries with the added benefit of being able to control major bleeding with the C-A-T® Tourniquet, the safest and most effective prehospital tourniquet on the market today.

FEATURES/SPECIFICATIONS

- Workplace first aid kits that exceed ANSI/ISEA Z308.1-2021 Standard for Class A kits
- Able to treat the most common workplace injuries with the added benefit of being able to control major bleeding with the safest and most effective Tourniquet on the market
- Designed so the kit can be easily carried to the patient's side

... With the added benefit of major bleeding control with the safest, most effective tourniquet on the market today!



DIMENSIONS

- L 14.75 in. x H 8 in. x D 4 in.
- Weight: 2.94 lbs.

TRAUMA & FIRST AID KITS - CLASS A

ITEM #	OPTION
80-0947	w/Compressed Gauze
80-1045	w/Bleeding Control Dressing
80-1477	w/CELOX™ gauze



Learn more about this and other Products with a Mission® at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277


TRAUMA & FIRST AID KITS (TFAK) - CLASS A




PLUS+
one of these:



- OR -



- OR -



**Exceeds ANSI/ISEA Z308.1-2021
Standard for Class A Kits**

KIT CONTENTS

- 1 x C-A-T® Tourniquet, Orange
- 1 x HyFin® Vent Chest Seal Twin Pack
- 1 x NAR® Compressed Gauze
- 1 x 4 in. NAR® Responder ETD®
- 1 x 4 in. Elastic Wrap Bandage
- 1 x 2 in. Elastic Wrap Bandage
- 1 x Triangular Bandage
- 1 x SAM® Splint
- 2 x Gauze, Pad - 4 in. x 4 in. (Pack of 2)
- 2 x Gauze Pad - 5 in. x 9 in.
- 16 x Flexible Fabric Bandages, 1 in. x 3 in.
- 1 x BurnTec® Dressing - 4 in. x 4 in.
- 1 x CPR Microshield Mask
- 10 x Antiseptic Towelette
- 10 x Antiseptic Ointment
- 1 x Roll 1 in. Surgical Tape
- 1 x Instant Cold Pack - 6 in. x 9 in.
- 1 x Eye Wash Solution, 1 oz
- 2 x Sterile Oval Eye Pad
- 2 pr Responder Nitrile Gloves - Large
- 1 pr Responder Trauma Shears - Large
- 10 x Hand Sanitizer
- 10 x Burn Gel
- 1 x NAR® Survival Blanket
- 1 x Black Permanent Marker, Large
- 1 x First Aid Pocket Guide
- 1 x Just In Time Bleeding Control Instructions
- 1 x Trauma & First Aid Kit Contents Inventory Card
- 1 x Inspection Card

Plus one of these:

- 1 x NAR® Compressed Gauze
- 1 x QuikClot® Bleeding Control Dressing
- 1 x Celox™ First Aid Hemostatic Dressing

Learn more about this and
other Products with a Mission®
at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277