

# INDIVIDUAL AID KIT

North American Rescue's compact, vacuum-sealed Individual Aid Kit is designed to treat bleeding from penetrating and other traumatic injuries. The items in this kit are simple to use by non-medically trained caregivers. The purpose of this kit is to control bleeding and cover penetrating injuries until medically trained responders arrive on scene.

**Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)**




## FEATURES/SPECIFICATIONS

- Compact first aid kit to treat penetrating and other traumatic injuries
- Designed for non-medically trained caregivers to control bleeding until medically trained responders arrive on scene
- Vacuum packaged & easy-to-open with the medical supplies needed to begin treatment

## DIMENSIONS

- H 5.75 in. x W 3.5 in. x D 2 in.
- Weight: 6.6 oz



## all KITS INCLUDE:



## KIT CONTENTS

- 1x NAR® Compressed Gauze
- 1x Responder Flat ETD, Emergency Trauma Dressing
- 1x Petrolatum Gauze
- 1x 1 in. Surgical Tape
- 1pr Responder Nitrile Gloves, Large

## IND. AID KIT

ITEM #

85-0404