

LOGGING TRAUMA & FIRST AID KIT (SOFT CASE)

"By many measures, logging is the most dangerous occupation in the United States."

"Logging operations are often remote from any emergency medical care. Properly administered first aid can be a major factor in minimizing even life-threatening injuries."

- OSHA.org

The Logging Trauma & First Aid Kit from North American Rescue is designed as a portable First Aid Kit for the workplace that exceeds OSHA Standard Number: 1910.266 App A - First-aid Kits (Mandatory) for the logging industry.

Packaged in a convenient zippered nylon bag that is easily transported to the patient's side, this kit offers the ability to treat logging environment injuries with the added benefit of being able to control major bleeding from traumatic injuries with the C-A-T® Tourniquet, the fastest, safest, most effective prehospital field tourniquet on the market today.

DIMENSIONS

- Packaged: L 14.75 in. x H 8 in. x D 4 in.
- Weight: 3.03 lbs

FEATURES/SPECIFICATIONS

- Compact size allows for storage in vehicles and machinery
- Rugged nylon constructed bag keeps contents secure
- Durable 550 cord pulls
- Able to treat the most common job site injuries with the added benefit of being able to control major bleeding with the fastest, safest, most effective prehospital field tourniquet on the market today
- Designed with mobility in mind, the kit can be easily carried to the patient's side



Learn more about this and other **Products with a Mission®** at www.NARescue.com



LOGGING TRAUMA & FIRST AID KIT (SOFT CASE)

ITEM #	OPTION
80-1188	Basic
80-1189	w/Bleeding Control Dressing
80-1910	w/ CELOX™




NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

LOGGING TRAUMA & FIRST AID KIT (SOFT CASE)

all KITS INCLUDE:



+plus,
by kit type:



KIT CONTENTS

- 1 x C-A-T® Tourniquet, Orange
- 1 x NAR® Responder ETD®, 6 in. Emergency Trauma Dressing
- 2 x Gauze, Pad, 8 in. x 10 in. (Pack of 2)
- 2 x Gauze, Pad, 4 in. x 4 in. (Pack of 2)
- 2 x Elastic Wrap Bandage, 2 in.
- 2 x Triangular Bandage
- 16 x Bandages, 1 in. x 3 in.
- 1 x BurnTec® Dressing, 4 in. x 4 in.
- 1 x SAM® Splint II
- 2 pr Responder Nitrile Gloves - Blue, Large
- 10 x Antiseptic Towelette
- 10 x Sting & Bite Towelette
- 1 x CPR Microshield Mask
- 1 x Eye Wash Solution, 1 oz
- 2 x Polycarbonate Eye Shield - Orange
- 2 x Instant Cold Compress
- 1 x Mini Surgical Tape, 2 in.
- 1 x Splinter Forceps, 4.5 in.
- 1 x Responder Shears, Large
- 1 x NAR® Survival Blanket
- 1 x Inspection Card
- 1 x Inventory/Content Card
- 1 x Just In Time Bleeding Control Instructions
- 1 x First Aid Pocket Guide

+plus, by kit type:

- 2 x NAR® Compressed Gauze -or-
- 1 x NAR® Compressed Gauze
- 1 x QuikClot® Z-Fold Bleeding Control Dressing -or-
- 1 x NAR® Compressed Gauze
- 1 x CELOX™ Z-Fold Hemostatic Dressing - 3 in. x 5 yd.



Learn more about this and other **Products with a Mission®** at www.NARescue.com



LOGGING TRAUMA & FIRST AID KIT (SOFT CASE)

ITEM #	OPTION
80-1188	Basic
80-1189	w/Bleeding Control Dressing
80-1910	w/ CELOX™

* Explore other Logging Trauma Kit Options at www.NARescue.com.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277