

LOGGING FIRST AID & TRAUMA KIT (HARD CASE)

"By many measures, logging is the most dangerous occupation in the United States."

"Logging operations are often remote from any emergency medical care. Properly administered first aid can be a major factor in minimizing even life-threatening injuries." - OSHA.org

The Logging First Aid & Trauma Kit from North American Rescue® is designed as a portable First Aid Kit for the workplace that exceeds OSHA Standard Number: 1910.266 App A - First-aid Kits (Mandatory) for the logging industry. Packaged in a durable dual-latched hard case that is easily transported to the patient's side, this kit offers the ability to treat logging environment injuries with the added benefit of being able to control major bleeding from traumatic injuries with the C-A-T® Tourniquet, the fastest, safest, most effective prehospital field tourniquet on the market today.

FEATURES/SPECIFICATIONS

- Rugged case designed for use in the harshest of conditions
- Modular system featuring a rugged storage organizer in the lid and removable nylon pouches with easy-to-read treatment labels (Meets Type IV ANSI Container)
- Fast open dual-latch closure system allows for rapid access while keeping contents secure
- Able to treat the most common job site injuries with the added benefit of being able to control major bleeding with the safest and most effective Tourniquet on the market
- Designed with mobility in mind, the kit can be easily carried to the patient's side

LOGGING TRAUMA & FIRST AID KITS - HARD CASE

ITEM #	COLOR	OPTION
80-1190	RED	Basic
80-1191	RED	w/Bleeding Control Dressing
80-1899	RED	w/CELOX™ gauze
80-1894	ORG	Basic
80-1896	ORG	w/Bleeding Control Dressing
80-1898	ORG	w/CELOX™ gauze
80-1895	YLW	Basic
80-1897	YLW	w/Bleeding Control Dressing
80-1900	YLW	w/CELOX™ gauze



Learn more about this and other Products with a Mission® at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

LOGGING FIRST AID & TRAUMA KIT (HARD CASE)



KIT CONTENTS

- 1 x C-A-T® Tourniquet, Orange
 - 1 x NAR® Responder ETD®, 4 in. Emergency Trauma Dressing
 - 2 x Gauze, Pad, 8 in. x 10 in. (Pack of 2)
 - 2 x Gauze, Pad, 4 in. x 4 in. (Pack of 2)
 - 2 x Elastic Wrap Bandage, 2 in.
 - 2 x Triangular Bandage
 - 16 x Bandages, 1 in. x 3 in.
 - 1 x BurnTec® Dressing, 4 in. x 4 in.
 - 1 x SAM® Splint II
 - 2 pr Responder Nitrile Gloves - Blue, Large
 - 10 x Antiseptic Towelette
 - 10 x Sting & Bite Towelette
 - 1 x CPR Microshield Mask
 - 1 x Eye Wash Solution, 1 oz
 - 2 x Polycarbonate Eye Shield - Orange
 - 1 x Instant Cold Compress
 - 1 x Mini Surgical Tape, 2 in.
 - 1 x Splinter Forceps, 4.5 in.
 - 1 x Responder Shears, Large
 - 1 x NAR® Survival Blanket
 - 1 x Inventory/Content Card
 - 1 x Just In Time Bleeding Control Instructions
 - 1 x First Aid Pocket Guide
- +plus, by kit type:**
- 2 x NAR® Compressed Gauze
 - -OR-
 - 1 x NAR® Compressed Gauze
 - 1 x QuikClot® Z-Fold Bleeding Control Dressing
 - -OR-
 - 1 x NAR® Responder Compressed Gauze
 - 1 x Celox™ First Aid Z-Fold Hemostatic Dressing - 3 in. x 5 ft.



Learn more about this and other **Products with a Mission®** at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277