

# AMPHIBIOUS AID ADVANCED KITS

## WATERTIGHT

### Keep Sensitive Gear Dry While Conducting Maritime And Underwater Operations

The Amphibious Aid Kit (Watertight)- from North American Rescue® is designed as a portable First Aid Kit for use in the workplace that exceeds the ANSI/ISEA Z308.1-2021 Standard for Class B kits. This is accomplished by adding additional bleeding control elements proven to save lives. To bring this kit to the market, NAR® partnered with Watershed®, the manufacturer of drybags and other watertight containers used by the US military special operations community to keep their sensitive gear dry while conducting maritime and underwater operations. The kit contents are arranged in labeled and zipper-accessed mesh pouches. This kit supplies not only the item needed to treat the most common workplace injuries but many of the most serious, life threatening with the C-A-T® Tourniquet, the safest and most effective prehospital tourniquet on the market today. This kit comes in Advanced with Combat Gauze®, and Advanced with Celox™. The Advanced versions additionally include a NPA and BVM for breathing capabilities.

**Learn more about this and  
other Products with a Mission®  
at [www.NARescue.com](http://www.NARescue.com)**



### FEATURES/SPECIFICATIONS

- A portable First Aid Kit for use in the workplace that exceeds the ANSI/ISEA Z308.1-2021 Standard for Class B kits.
- NAR® partnered with Watershed to keep their sensitive gear dry while conducting maritime and underwater operations.
- Supplies not only the item needed to treat the most common workplace injuries but many of the most serious, life threatening ones
- Comes in Advanced with Combat Gauze® and Advanced with Celox™. The Advanced versions additionally include a NPA and BVM for breathing capabilities.



### AMPHIBIOUS AID ADV. KITS - (WATERTIGHT)

ITEM #	COLOR	OPTION	NSN
80-1617	Orange	Advanced with Combat Gauze®	----
80-1618	Ranger Green	Advanced with Combat Gauze®	----
80-1619	Orange	Advanced with Celox™ Gauze	----
80-1620	Ranger Green	Advanced with Celox™ Gauze	----

\* Explore basic versions of this kit at [www.NARescue.com](http://www.NARescue.com).

### DIMENSIONS

- H 6.5 in. x W 9 in. x L 11.5 in.
- Weight: 4 lbs 11.2 oz



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

# AMPHIBIOUS AID ADVANCED KITS

## WATERTIGHT

### all KITS INCLUDE:



### +plus, by kit type:



## KIT CONTENTS

- 2 x C-A-T® - Orange (Combat Application Tourniquet®)
- 2 x ETD® 4 in Flat Responder Emergency Trauma Dressing
- 1 x HyFin® Vent Compact Chest Seal - Twin Pack
- 2 x NAR® Wound Packing Gauze - Z-Fold 3 in x 5 yd Dressing
- 1 x Pad Gauze - 4 in x 4 in (2 Pack)
- 2 x Pad Gauze - 5 in x 9 in
- 2 x NAR® Pre-lubricated Nasopharyngeal Airway - 28FR
- 1 x Bag Valve Mask (BVM) - Cyclone
- 1 x Instant Cold Pack - 6 in x 9 in
- 1 x Universal Malleable Splint
- 1 x NAR® Survival Blanket - 56 in. x 96 in.
- 5 pr Sand Bear Claw® Nitrile Glove Kit
- 1 x 7.25 in. Responder Trauma Shears
- 2 x Sterile Oval Eye Pad - 1.625 in x 2.625 in
- 1 x Gecko Grip Multi-purpose Flat Tape (6 Count Strips - 2" x 15")
- 2 x Celox™ Gauze - Z-fold 3 in x 5ft Dressing
- 2 x Elastic Wrap Bandage - 4 in 5 yd
- 1 x Triangular Bandage - 40 in x 40 in x 56 in
- 1 x Elastic Adhesive Tape - 4 in 5 yd
- 1 x Large Permanent Marker, Black
- 1 x MED Cross Patch (2"x2") - Black/Green

Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



## AMPHIBIOUS AID ADV. KITS - (WATERTIGHT)

ITEM #	COLOR	OPTION	NSN
80-1617	Orange	Advanced with Combat Gauze®	----
80-1618	Ranger Green	Advanced with Combat Gauze®	----
80-1619	Orange	Advanced with Celox™ Gauze	----
80-1620	Ranger Green	Advanced with Celox™ Gauze	----

\* Explore basic versions of this kit at [www.NARescue.com](http://www.NARescue.com).



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277