

# SMRT™ TASK TRAINER - PACKABLE FEMORAL WOUND WITH GUNSHOT/BLAST

## STRATEGIC MEDICAL RESEARCH & TRAINING™

The Packable Femoral Wound with Gunshot / Blast wearable injury from SMRT™ empowers users to experience treatment of anatomically accurate injury patterns during instructional presentations, live action scenarios, and trauma lane practice sessions. This injury features a bleeding femoral wound that responds to direct pressure and wound packing—allowing students to experience instant treatment confirmation during practice.

Large cavity, high femoral wounds are among the most challenging, life-threatening wounds. This wound simulates a blast or gunshot exit wound capable of accepting five rolls of gauze. The unique design prevents slipping while worn and anatomically corrects bleeding.

The innovative and wearable SMRT™ Moulage Task Trainers from Strategic Medical Research & Training™ create an unmatched experience of realism featuring realistic moulage and active bleeding across a wide spectrum of woundings and interventions. Made with durable materials and designed for a "one size fits most" range of limb and body sizes, the SMRT™ Moulage Task Trainers are engineered to work with live role players as well as to augment existing simulation systems—while allowing flexibility in placement and utilization. All SMRT™ Moulage Task Trainers work with SMRT™ Trauma Tables, SMRT™ Bleeding Systems, and common hydration systems, making them standalone and ready to use, while also allowing ease of integration into existing training systems. Utilized by Military, Law Enforcement, EMS, and Hospitals to create unsurpassed realism in hemorrhage control training, SMRT™ Training Aids are realistic, cost effective, highly mobile, durable, and low maintenance.



**Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)**



## FEATURES/SPECIFICATIONS

- Unsurpassed realism in hemorrhage control training
- One size fits most range of limb and body sizes
- Packable Femoral Wound with Gunshot/Blast compatible with other SMRT™ Task Trainers
- Available in Light Tone or Deep Tone

## SMRT™ PACKABLE FEMORAL - GSW/BLAST

ITEM #	COLOR
93-1183	Light Tone
93-1184	Deep Tone

\* Explore other SMRT™ moulage task trainers at [www.NARescue.com](http://www.NARescue.com)

## DIMENSIONS

- Packaged Dimensions: L 9 in. x W 11 in. x D 3 in.
- Weight: 2.6 lbs.



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277