

# NANO IFAK

The Nano IFAK is the result of a collaboration between North American Rescue® and Blue Force Gear®, featuring NAR's premium trauma first aid components with BFG's cutting edge custom designed bag. The Nano IFAK features rapid single-hand access by simply pulling the BLIP pull tab to remove the insert and access the medical supplies. Additionally, the NANO bag is designed with an external tourniquet holder and Naloxone pouch (Naloxone not included), allowing swift access to critical life saving gear without opening the main compartment. With the ability to mount vertically or horizontally on a belt or MOLLE platform, it delivers the perfect marriage of optimized cube space with life saving trauma kit capabilities.

The Nano IFAK has two mounting options – featuring a belt mounting slot for a duty belt or MOLLE for compatible plate carriers or other gear:

- The Belt mount features a simple, lightweight, and durable way to quickly attach the Nano IFAK to any duty belt or pants belt 2 in. or smaller. Constructed out of Blue Force Gear's ULTRAcamp high-performance coated fabric laminate, the backer offers a stiff yet low-profile and reliable mounting solution for either horizontal or vertical carry.
- The MOLLE attachment features Blue Force Gear's two-way Helium Whisper attachment system that allows for pouch mounting both vertically and horizontally. BFG's patented Helium Whisper attachment system is the lightest and strongest MOLLE style attachment system available and is the choice of the most elite professionals.

**Learn more about this and  
other Products with a Mission®  
at [www.NARescue.com](http://www.NARescue.com)**



## FEATURES/SPECIFICATIONS

- A collaboration with North American Rescue® and Blue Force Gear®
- A compact, lightweight trauma kit with full size capabilities
- Kit contains the ability to mount vertically or horizontally on a belt or MOLLE platform
- Features rapid single-hand access by simply pulling the BLIP pull tab to remove the insert and access medical supplies
- Nylon bag is designed with built-in C-A-T® Tourniquet holder and Naloxone pouch (Naloxone not included)



**Blue Force Gear bag  
features Naloxone  
Pouch (Naloxone  
not included)**

## NANO IFAK

ITEM #	COLOR	OPTION	ATT.
80-1942	Black	Basic	Belt
80-1943	Ranger	Basic	Belt
80-1946	Black	Basic	MOLLE
80-1947	Ranger	Basic	MOLLE
80-1944	Black	w/ Combat Gauze®	Belt
80-1945	Ranger	w/ Combat Gauze®	Belt
80-1948	Black	w/ Combat Gauze®	MOLLE
80-1949	Ranger	w/ Combat Gauze®	MOLLE

\* Similar kits containing **Naloxone** are available at [www.NARescue.com](http://www.NARescue.com).

## DIMENSIONS

- Kit Dimensions:
  - » Full Kit (with C-A-T® Tourniquet):  
H 7 in. x W 5.5 in. x D 2 in.
  - » Not including C-A-T®:  
H 5 in. x W 5.5 in. x D 2 in.
  - » With Kit in resealable bag:  
H 11 in. x W 9 in. x D 2 in.
- Weight: 0.82 lb.



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

# NANO IFAK

## all KITS INCLUDE:



## +plus, by kit type:



## KIT CONTENTS

- 1x Combat Application Tourniquet® (C-A-T®) - Black
- 1x NAR® Mini Responder ETD®, (Emergency Trauma Dressing) - 4 in.
- 1x Hyfin® Vent Compact Chest Seal (Twin Pack)
- 1pr Blue Responder Gloves - Large

## +plus, by kit type:

- 1x NAR® Wound Packing Gauze - 3 in. x 5 yd.
- or -
- 1x Combat Gauze® LE - 3 in. x 4 yd.



▲  
**BELT  
MOLLE**  
▼

## NANO IFAK

ITEM #	COLOR	OPTION	ATT.
80-0460	Black	Basic	Belt
80-0461	Ranger	Basic	Belt
80-0460	Black	Basic	MOLLE
80-0461	Ranger	Basic	MOLLE
80-0460	Black	w/ Combat Gauze®	Belt
80-0461	Ranger	w/ Combat Gauze®	Belt
80-0460	Black	w/ Combat Gauze®	MOLLE
80-0461	Ranger	w/ Combat Gauze®	MOLLE

\* Similar kits containing **Naloxone** are available at [www.NARescue.com](http://www.NARescue.com).



Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277