

# NANO IFAK - WITH NALOXONE

The Nano IFAK is the result of a collaboration between North American Rescue® and Blue Force Gear®, featuring NAR's premium trauma first aid components with BFG's cutting edge custom designed bag. The Nano IFAK features rapid single-hand access by simply pulling the BLIP pull tab to remove the insert and access the medical supplies. Additionally, the NANO bag is designed with an external tourniquet holder and Naloxone pouch, allowing swift access to critical life saving gear without opening the main compartment. Additionally, this version of the Nano IFAK comes with Naloxone™ as a pack of two and in its original packaging. With the ability to mount vertically or horizontally on a belt or MOLLE platform, it delivers the perfect marriage of optimized cube space with life saving trauma kit capabilities.

The Nano IFAK has two mounting options – featuring a belt mounting slot for a duty belt or MOLLE for compatible plate carriers or other gear.

- The Belt mount features a simple, lightweight, and durable way to quickly attach the Nano IFAK to any duty belt or pants belt 2 in. or smaller. Constructed out of Blue Force Gear's ULTRAcamp high-performance coated fabric laminate, the backer offers a stiff yet low-profile and reliable mounting solution for either horizontal or vertical carry.
- The MOLLE attachment features Blue Force Gear's two-way Helium Whisper attachment system that allows for pouch mounting both vertically and horizontally. BFG's patented Helium Whisper attachment system is the lightest and strongest MOLLE style attachment system available and is the choice of the most elite professionals.



**Blue Force Gear bag features Naloxone Pouch (Naloxone comes in a pack of 2 and in its original packaging)**

**Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)**



## FEATURES/SPECIFICATIONS

- A collaboration with North American Rescue® and Blue Force Gear®
- A compact, lightweight trauma kit with full size capabilities
- Kit contains the ability to mount vertically or horizontally on a belt or MOLLE platform
- Features rapid single-hand access by simply pulling the BLIP pull tab to remove the insert and access medical supplies
- Nylon bag is designed with built-in C-A-T® Tourniquet holder and Naloxone pouch (Naloxone comes in a pack of two, and in its original packaging)

## NANO IFAK - WITH NALOXONE

ITEM #	COLOR	OPTION	ATT.
80-1952	Black	Basic	Belt
80-1953	Ranger	Basic	Belt
80-1956	Black	Basic	MOLLE
80-1957	Ranger	Basic	MOLLE
80-1954	Black	w/ Combat Gauze®	Belt
80-1955	Ranger	w/ Combat Gauze®	Belt
80-1958	Black	w/ Combat Gauze®	MOLLE
80-1959	Ranger	w/ Combat Gauze®	MOLLE

\* Similar kits **without Naloxone** are available at [www.NARescue.com](http://www.NARescue.com).

## DIMENSIONS

- Kit Dimensions:
  - » Full Kit (with C-A-T® Tourniquet): H 7 in. x W 5.5 in. x D 2 in.
  - » Not including C-A-T®: H 5 in. x W 5.5 in. x D 2 in.
  - » With Kit in resealable bag: H 11 in. x W 9 in. x D 2 in.
- Weight: 0.91 lb.



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

# NANO IFAK - WITH NALOXONE

## all KITS INCLUDE:



## +plus, by kit type:



## KIT CONTENTS

- 1x Combat Application Tourniquet® (C-A-T®) - Black
- 1x NAR® Mini Responder ETD®, (Emergency Trauma Dressing) - 4 in.
- 1x Hyfin® Vent Compact Chest Seal (Twin Pack)
- 1pr Blue Responder Gloves - Large
- 1x Naloxone™ (2pk)

## +plus, by kit type:

- 1x NAR® Wound Packing Gauze - 3 in. x 5 yd.
- or -
- 1x Combat Gauze® LE - 3 in. x 4 yd.



▲  
**BELT  
MOLLE**  
▼

## NANO IFAK - WITH NALOXONE

ITEM #	COLOR	OPTION	ATT.
80-1952	Black	Basic	Belt
80-1953	Ranger	Basic	Belt
80-1956	Black	Basic	MOLLE
80-1957	Ranger	Basic	MOLLE
80-1954	Black	w/ Combat Gauze®	Belt
80-1955	Ranger	w/ Combat Gauze®	Belt
80-1958	Black	w/ Combat Gauze®	MOLLE
80-1959	Ranger	w/ Combat Gauze®	MOLLE

\* Similar kits **without** Naloxone are available at [www.NARescue.com](http://www.NARescue.com).



Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277