

OUT-PAK (LCL)

Outdoor Personal Aid Kit (Laser Cut Laminate)

A compact Mini First Aid Kit...now in LCL Nylon!

North American Rescue's Outdoor Personal Aid Kit (OUT-PAK) - now featuring Laser Cut Laminate (LCL) MOLLE - was designed as a multi-use IFAK for delivering immediate care for injuries encountered in an outdoor adventure. At half the size of our popular TORK platform, the OUT-PAK is currently one of the smallest platforms available in our product offering. Despite its reduced size, this kit comes fully loaded with the essential medical equipment for Basic, Basic with Celox™ Hemostatic Dressing, or Basic with QuikClot® Bleeding Control Dressing.

This Individual First Aid Kit's super compact, rugged nylon bag opens via 550 cord zipper pulls from the side in a clam shell configuration utilizing (2) two main sleeves that open on both ends for easy access and features multiple elastic loops for secure gear organization. The bag's exterior includes a MOLLE backing for mounting in the standard vertical position as well as a 3 inch internal sleeve for accommodating horizontal attachment to your duty belt. The vertical mount can be set to open left-to-right or right-to-left based on shooter preference, while the horizontal mount on a belt allows opening directly to your C-A-T® tourniquet. Also included are a loop patch on front of the bag for custom user labeling and (3) three MOLLE strips on the front, and (2) two MOLLE strips on the spine's edge to attach additional gear such as TQ/holders, tactical lights, or knives (sold separately) if needed.

FEATURES/SPECIFICATIONS

- NAR's most compact, versatile Outdoor Personal Aid Kit Platform designed for outdoor adventure
- Contains essential point-of-wounding medical equipment for treating injuries
- Half the size of our popular TORK platform
- The LCL IFAKS are double spaced for the MOLLE attachment points on the rear of each IFAK for maximum adjustability and super secure mounting to any platform
- Compatible with NAR® Quick Detach Panels (sold separately)

Learn more about this and other Products with a Mission® at www.NARescue.com



OUT-PAK KIT

COLOR	LEVEL	ITEM#	NSN #
BLK	Basic	80-1821	----
COY	Basic	80-1822	----
OCP	Basic	80-1823	----
ODG	Basic	80-1824	----
RGR	Basic	80-1826	----
RED	Basic	80-1825	----
BLK	w/Combat Gauze®	80-1833	----
COY	w/Combat Gauze®	80-1834	----
OCP	w/Combat Gauze®	80-1835	----
ODG	w/Combat Gauze®	80-1836	----
RGR	w/Combat Gauze®	80-1838	----
RED	w/Combat Gauze®	80-1837	----
BLK	w/CELOX™	80-1827	----
COY	w/CELOX™	80-1828	----
OCP	w/CELOX™	80-1829	----
ODG	w/CELOX™	80-1830	----
RGR	w/CELOX™	80-1832	----
RED	w/CELOX™	80-1831	----

DIMENSIONS


- L 6 in. x W 3 in. x D 3 in.
- Weight: 0.99 lbs.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

OUT-PAK (LCL)

Outdoor Personal Aid Kit (Laser Cut Laminate)



KIT CONTENTS

- 1 x C-A-T® (Combat Application Tourniquet®) - ORG
- 1 x NAR® Mini Responder ETD®, 4in. (Emergency Trauma Dressing)
- 1 x NAR® Survival Blanket
- 5 x Flexible Fabric Bandage (3/4 in. x 3 in.)
- 1 x Wound Closure Strip Bandage (pack of 10)
- 3 x Antibiotic Ointment
- 2 x Sting & Bite Packet
- 3 x Safetec .9 G Burn Gel
- 1 x Mini Duct Tape (2 in. x 100 in.)
- 1 pr Bear Claw® Nitrile Gloves - SND, Lg
- 1 x Med Cross Patch - GRY/BLK/TAN, 2 in. x 2 in.

Plus one of these:

- 1 x NAR® Compressed Gauze
- 1 x QuikClot® Bleeding Control Dressing - 3 in. x 4 yd.
- 1 x CELOX™, Z-Fold, Hemostatic Gauze



Learn more about this and other **Products with a Mission®** at www.NARescue.com



OUT-PAK KIT

COLOR	LEVEL	ITEM#	NSN #
BLK	Basic	80-1821	----
COY	Basic	80-1822	----
OCP	Basic	80-1823	----
ODG	Basic	80-1824	----
RGR	Basic	80-1826	----
RED	Basic	80-1825	----
BLK	w/Combat Gauze®	80-1833	----
COY	w/Combat Gauze®	80-1834	----
OCP	w/Combat Gauze®	80-1835	----
ODG	w/Combat Gauze®	80-1836	----
RGR	w/Combat Gauze®	80-1838	----
RED	w/Combat Gauze®	80-1837	----
BLK	w/CELOX™	80-1827	----
COY	w/CELOX™	80-1828	----
OCP	w/CELOX™	80-1829	----
ODG	w/CELOX™	80-1830	----
RGR	w/CELOX™	80-1832	----
RED	w/CELOX™	80-1831	----



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277