

TRAUMA AND FIRST AID BOATING KIT

When major bleeding happens, seconds can be the difference!

Meet and exceed the United States Coast Guard (USCG) requirements for first aid kits on recreational boats. Don't be caught on the water without the proper medical kit in your boat, and don't be caught in an emergency without advanced life saving bleeding control solutions.

The Trauma and First Aid Boating Kit from North American Rescue® was built to exceed the recreational boating standards set by the United States Coast Guard. This kit includes the tools to care for small cuts, burns, eye injuries as well as critical life saving equipment to render emergency care for major hemorrhages while out on the water. Treat aches and other minor injuries with the included over-the-counter medications. The First Aid Pocket Guide walks care givers step-by-step through rendering first aid care to the injured party.

Kit items are packed into zippered nylon storage bags that offer the ability to organize injury related products for fast deployment. The red bag labels eliminates confusion when you need to quickly locate the proper treatment supplies. With a brightly colored exterior, the rugged water-tight case is easy to locate and keeps all component bags secure. The integrated carry handle provides a comfortable grip and enables simple movement to get the kit where it's needed for use.

FEATURES/SPECIFICATIONS

- Exceeds the U.S. Coast Guard Standard 160.041 - Items for First Aid Kits
- Treat a range of injuries from minor to massive hemorrhage
- Includes over-the-counter (OTC) medications
- Airtight, Waterproof and Sealable
- Rugged, easy to see hard case
- The pouches are designed to easily and quickly access the products necessary to respond to emergencies while on the water
- Contains one utility organizer mounted in the lid of the case



Learn more about this and other Products with a Mission® at www.NARescue.com



TRAUMA AND FIRST AID BOATING KIT

ITEM #	COLOR	OPTION
80-1886	Orange	Basic
80-1887	Red	Basic
80-1483	Yellow	Basic
80-1888	Orange	w/ Bleeding Control Dressing
80-1889	Red	w/ Bleeding Control Dressing
80-1890	Yellow	w/ Bleeding Control Dressing
80-1891	Orange	w/ CELOX™
80-1892	Red	w/ CELOX™
80-1893	Yellow	w/ CELOX™

DIMENSIONS


- Packaged: H 6 in. x W 14 in. x D 10 in.
- Weight: 9.4 lbs.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

TRAUMA AND FIRST AID BOATING KIT

all KITS INCLUDE:



Learn more about this and other Products with a Mission® at www.NARescue.com



KIT CONTENTS

- 1x C-A-T® Tourniquet - Orange
 - 1x NAR® Compressed Gauze
 - 1x NAR® Responder Emergency Trauma Dressing (ETD®) - 6 in.
 - 1x CPR Pocket Resuscitator
 - 1x HyFin® Vent Chest Seal Twin Pack
 - 4x Petrolatum Gauze - 3 in. x 18 in.
 - 1x Emergency Survival Wrap - L 96 in. x W 60 in.
 - 1x BurnTec Dressing - 4 in. x 4 in.
 - 5x Safetec Burn Gel
 - 3x Gauze Pad (2pk) - 4 in. x 4 in.
 - 1x Elastic Bandage - 4 in.
 - 4x Elastic Bandage - 2 in.
 - 1x SAM® Splint II
 - 2x Triangular Bandage - 37 in. x 37 in. x 52 in.
 - 1x Eye Wash Solution - 4 oz.
 - 1x Surgical Tape - 1 in.
 - 3x PES Eye Shield
 - 4x Saljet Rinse - 1 oz.
 - 1x Wound Closure Strip (10pk)
 - 10x Fabric Bandage - 1 in. x 3 in.
 - 5x Flexible Bandage - 2 in. x 4 in.
 - 10x Povidone-Iodine Anti-septic Prep Pad
 - 5x Antibiotic Ointment
 - 1x Bacitracin Ointment
 - 2x Antacid + Acid Reducer Chewable Tablet
 - 2x Over-the-counter (OTC) Ibuprofen
 - 2x OTC Diphenhydramine
 - 5x OTC Aspirin
 - 10x Ammonia Inhalant
 - 1x Splinter Forceps - 4.5 in.
 - 1x Curved Kelly Hemostat - 5.5 in.
 - 12x Safety Pins - 1 in.
 - 4pr Responder Blue Gloves - Large
 - 1x Responder Shears - 7.25 in.
 - 1x Dikes/Pliers - 6 in.
 - 1x First Aid Pocket Guide
 - 1x Permanent Marker - Large
- +plus, by kit type:**
- 1x NAR® Compressed Gauze
 - or -
 - 1x QuikClot® Bleeding Control Dressing (Z-Fold) - 3 in. x 4 yd.
 - or -
 - 1x Celox™ First Aid Z-Fold Hemostatic Dressing - 3 in. x 5 ft.

TRAUMA AND FIRST AID BOATING KIT

ITEM #	COLOR	OPTION
80-1886	Orange	Basic
80-1887	Red	Basic
80-1483	Yellow	Basic
80-1888	Orange	w/ Bleeding Control Dressing
80-1889	Red	w/ Bleeding Control Dressing
80-1890	Yellow	w/ Bleeding Control Dressing
80-1891	Orange	w/ CELOX™
80-1892	Red	w/ CELOX™
80-1893	Yellow	w/ CELOX™



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277