

# LEATHERMAN RAPTOR EMERGENCY MEDICAL MULTITOOL

Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



## Optional Action Statement or Major Fact

North American Rescue® is proud to offer the latest innovation from Leatherman®, the industry leader in application-driven, quality-made multi-tools for the last 30 years. With the right balance of features for medical-specific ops without an overkill of options, the Leatherman® Raptor™ is a high-performance multi-tool equipped with a variety of mission critical tools essential to emergency medical practitioners in the field. Components include precision 420HC stainless steel medical shears with a stow-away strap cutter and built-in oxygen tank wrench that may be quickly accessed with one hand, and can be used when the shears are open or closed. Additional components include an ultra-strong ring cutter that will easily cut through most standard gold bands, a carbide-tip glass breaker and handy 5 cm ruler. The foldable shears come apart for easy cleaning and feature serrated blades to help grip materials for a strong fast cut. The ergonomic handles also feature glass-filled nylon scales and sport a lanyard hole and replaceable pocket clip for a variety of useful attachment points. The Raptor™ comes complete with an injection molded polymer holster which is also MOLLE-compatible and allows device storage in the open or folded position. The holster's secure rotating clip easily attaches to your vest, pack or duty belt allowing rapid access to your multi-tool when you need it. Backed by a 25 year warranty, the Leatherman® Raptor™ is built to last.. a must have for every mobile medic and 1st responder working in life-threatening emergency situations.



## LEATHERMAN RAPTOR EMERGENCY MEDICAL MULTITOOL

ITEM #	COLOR	NSN
70-0014	Black	5110-01-626-8766
70-1560	Blue/Black	----
70-1561	Blue/Yellow	----
70-1562	Coyote/Black	----



## FEATURES/SPECIFICATIONS

- Integrated multi-tool specific to emergency medical situations
- High quality 420HC stainless steel medical shears
- One-handed Strap Cutter
- Stow-away Oxygen Tank Wrench
- Injection-molded polymer Holster on a rotating clip; MOLLE-compatible

## DIMENSIONS

- Shear length: 1.9 in.
- Closed length: 5 in.
- Weight: 5.8 oz



# LEATHERMAN RAPTOR EMERGENCY MEDICAL MULTITOOL



1. 420HC Stainless Steel Folding Medical Shears
2. Strap Cutter
3. Oxygen Tank Wrench
4. Ring Cutter
5. Carbide Glass Breaker
6. 5 cm Ruler
7. Holster with Clip

Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)

