

TECHLINE TRAUMA BLOOD SUPPLY PUMP

The TechLine Trauma® Blood Supply Pump is a simulated blood supply reservoir with manual pump, designed to attach to any bleeding wearable wound simulators, Wound in a Box™ task trainers, or part task trainers.

The blood flow is controlled by the operator by squeezing the hand-bulb attached to the supply line. This allows the operator to stop the blood flow when the proper medical intervention has been taken.

CONSTRUCTION

- 2 liter hydration reservoir bag with screw top opening.
- Hand-bulb pump with check valve to prevent open flow.
- 30 inches of blood flow tubing with standard female connector.

**Learn more about this and
other Products with a Mission®
at www.NARescue.com**



FEATURES/SPECIFICATIONS

- Compatible with any bleeding wearable wound simulators, TechLine Trauma® Wounds in a Box™ task trainers, or part task trainers
- Manual pump allows operator to stop blood flow when proper medical intervention has been taken

BLOOD SUPPLY PUMP

ITEM #	NSN
93-1024	6910-01-672-4172

* Explore more TechLine Trauma products at www.NARescue.com

DIMENSIONS

- Packaged Measurements: L 15 in. x W 7 in.
- Weight: 0.63 lbs. lbs.