

TECHLINE TRAUMA TOMMANIKIN LEFT ARM

UNINJURED

The TechLine Trauma® TOMManikin® uninjured left arm features radial pulse, antecubital IV access (with flash), and humeral IO access (procedural only). The simulator's pulse rate, pressure, and regularity controlled via NARS tablet.

Features 100% silicone construction with metal bone internals. Interchangeable with all TechLine Trauma® TOMManikin® and TAMlkin® simulators.

Manufactured by Techline Trauma®

Learn more about this and
other *Products with a Mission*®
at www.NARescue.com



FEATURES/SPECIFICATIONS

- Features radial pulse
- Antecubital IV access
- Humeral IO access (procedural only)
- 100% silicone construction

TOMM LEFT ARM - UNINJURED

ITEM #	NSN
93-0088	----

* Explore more TechLine Trauma products at www.NARescue.com

DIMENSIONS

- L 24 in. x W 12 in. x H 6 in.
- Weight: 10 lbs.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277