

TECHLINE TRAUMA TOMMANIKIN GUNSHOT WOUND CHEST SKIN

The TechLine Trauma® TOMManikin® features anterior entry GSW with non-linear (post-axillary) exit wound. Actively bleeds from both wound locations. Anterior entry wound mixes with air from onboard compressor to create bubbling bleeding effect. Supports chest seal application and pleural needle decompression. Bilateral chest tube compatible (no additional charge, requires request at time of order).

Manufactured by Techline Trauma®

Learn more about this and
other Products with a Mission®
at www.NARescue.com



FEATURES/SPECIFICATIONS

- 100% silicone construction
- Interchangeable with all TechLine Trauma® TOMManikin® simulators
- Supports Chest Seal and Pleural Needle decompression
- Non-linear (post-axillary) exit wound

TOMMANIKIN® CHEST SKIN - GSW

ITEM #	NSN
93-0149	----

* Explore more TechLine Trauma products at www.NARescue.com

DIMENSIONS

- Packaged - L 20 in. x W 6 in. x H 12 in.
- Weight: 8 lbs



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277