

TECHLINE TRAUMA SIMULATOR IO HUMERAL BALL REPLACEMENT

The TechLine Trauma® TOMManikin® Humeral Ball Replacement allows for multiple training insertions. Humeral IO site is for procedural insertion only. This part does not support infusion. The Humeral Ball Replacement is interchangeable with all TechLine Trauma® TOMManikin® and TAMlkin® simulators.

Manufactured by Techline Trauma®

Learn more about this and other Products with a Mission® at www.NARescue.com



IO HUMERAL BALL REPL.

ITEM #	NSN
93-0178	----

* Explore more TechLine Trauma products at www.NARescue.com

FEATURES/SPECIFICATIONS

- Interchangeable with all TechLine Trauma® TOMManikin® / TAMlkin® simulators

DIMENSIONS

- Packaged - L 6 in. x W 6 in. x H 6 in.
- Weight: 2 lbs