

TECHLINE TRAUMA BLEEDING INTERVENTION TRAINER

SHRAPNEL

The TechLine Trauma® Shrapnel Bleeding Intervention Trainer is a mobile solution specifically designed for point-of-injury trauma simulation and hemorrhage control. The system incorporates a versatile fluid delivery mechanism, featuring a fluid bag that supports both gravity-fed and manual pump operation, allowing precise control over the delivery of simulated blood.

The TechLine Trauma® Shrapnel Wound in a Box® task trainer simulates a gunshot to soft-tissue, used to instruct students in proper wound packing technique. The representative skin and tissue of this product is constructed of high quality silicone rubber, with layers of coloration that provide a depth and realism to the product. The PATENTED construction with SOFTTECH® materials a realistic feel, and is very durable and UV resistant. Blood tubing is integrated into the construction to easily create simulated blood flow. The task trainer is contained in its own storage case.

Manufactured by TechLine Trauma®



FEATURES/SPECIFICATIONS

- Easily transportable, yet rugged enough for classroom and field use
- Set up/ Break down in less than a minute
- Highly Adaptive Bleeding Control
- Focused Hemorrhage Control Training
- Patented materials give this product a realistic feel, while still being durable and UV resistant

BLEEDING INTERVENTION TRAINER - SHRAPNEL

ITEM #	NSN
93-1118	----

* There are other options available at a place for this item and other items.

DIMENSIONS

- Measurements (Exterior): L 18.5 in. x W 14.06 in. x D 6.93 in. (47 cm. x 35.7 cm. x 17.6 cm.)
- Weight: 14 lbs.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

TECHLINE TRAUMA BLEEDING INTERVENTION TRAINER

SHRAPNEL

KIT INCLUDES:



CONSTRUCTION

The Shrapnel Wound in a Box® is contained in its own storage case. The task trainer reinforces wound packing technique for this unique mechanism of injury. Built within a durable Pelican 1500 case, this durable trainer offers two stands for fluid reservoirs to simulate massive hemorrhage.

OPERATION

The Shrapnel Wound in a Box® is contained in its own storage case. The task trainer reinforces wound packing technique for this unique mechanism of injury. Built within a durable Pelican 1500 case, this durable trainer offers two stands for fluid reservoirs to simulate massive hemorrhage.

MAINTENANCE

The device should be rinsed in clean, warm water and allowed to dry before storage. The simulation may be removed from the storage case for cleaning. The device may be washed with mild detergents like dish soap. Allow the device to dry fully before storage to prevent mold or mildew growth.

KIT CONTENTS

- 1 x Wound in a Box - Shrapnel
- 1 x Blood Supply Pump with Carry Bag
- 1 x Bladder Mounting Rod
- 2 x Simulated Blood - 8 oz.
- 8 x NAR® Wound Packing Gauze
- 2 x QuikClot® Combat Gauze Trainer
- 1 x Portable Bleed Tray
- 1 x Pelican Hard Case

BLEEDING INTERVENTION TRAINER - SHRAPNEL

ITEM #	NSN
93-1118	----

* There are other options available at a place for this item and other items.

Learn more about this and other Products with a Mission® at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277