

TCCC MAN[®] SIMULATOR

Learn more about this and other **Products with a Mission[®]** at www.NARescue.com

Innovative, Cutting-Edge Medical Training

The TOMManikin[®] / TAMiKin[™] (Tactical Operations Medical Manikin), co-developed by NAR Simulation[®] and Techline Technologies Inc., is an advanced training system which elevates Tactical Combat Casualty Care (TCCC) training from the point of injury through transfer to higher levels of care, delivering unmatched durability and versatility.

Includes GSW chest with anterior entry and mid-axillary exit wound, deep packable pelvic wound, complete below-the-knee amputation leg (TQ responsive), and unstable pelvis with palpatus crepitus



Realistic, Mission-Ready Simulation

The TOMM Series enables instructors to build highly realistic, traumatic casualty scenarios with responsive, reactive patients. Its robust design supports use in classrooms, trauma lanes, and demanding field training environments.

FEATURES / SPECIFICATIONS

- Supports realistic mission profiles for combat operations, disaster simulations, and mass casualty training

Durable, Rugged Construction

- Built with silicone and a welded steel frame
- Weighs **150 lbs.**
- Extremities withstand up to **6,000 lbs.** of pressure — ideal for confined spaces, collapsed structures, and care under fire

Advanced Remote Control & Feedback

- Includes **NAR Simulation Advanced Remote Tablet & Software** to control up to 12 units wirelessly
- Manages bleeding locations and bleeding type (venous or arterial)

Comprehensive Intervention Capabilities

- Tourniquet application (proximal, junctional, abdominal)
- Wound packing with bleeding control
- Airway management (NPA, OPA, supraglottic, intubation)
- Sternal and Tibial IO access
- Realistic chest seal placement
- CPR chest compressions

Field-Ready Power & Fluids

- Onboard 12V rechargeable NiMH batteries with chargers
- **1.2 gallon** (4.5 liters) onboard blood reservoir for multiple bleeding sites

3-year WARRANTY



Water-Resistant & Low Maintenance

- Designed for tough environments with quick setup, rapid reset, and minimal upkeep

DIMENSIONS

- Carrier Bags (2 x): L 42" x W 24" x H 16"
- Assembled Manikin: L 74" x W 24" x H 10"
- Weight: 150 lbs.

TCCC MAN[®]

ITEM #	NSN
93-0958	----

* Explore more simulator options at www.NARescue.com



NORTH AMERICAN RESCUE[®]
www.NARescue.com • 888.689.6277

TCCC MAN[®] SIMULATOR

Train Real. Train Ready.

From urban tactical operations to disaster response and combat casualty care, the TechLine Trauma[®] TOMManikin[®] / TAMikin[™] provides the realism, durability, and control you need to build confident, capable responders.

3-year WARRANTY

Learn more about this and other Products with a Mission[®] at www.NARescue.com



TCCC MAN[®] PACKAGE INCLUDES

- NAR Simulation Advanced Remote Tablet & Software
- 1.2 gallon onboard blood reservoir
- Pelican 1500 accessories case
- 12V NiMH batteries (10AH) with smart chargers (2x)
- T-handle tool & fluid connectors
- Hook & Loop training uniform (blouse, trousers, and belt)
- NAR[®] IFAK
- TOMManikin[®] Operator's Manual
- 3-wheeled deployment-style carrier bags for safe storage and transport (2x)



OPERATIONAL MEDICAL PROCEDURES

- Massive hemorrhage control (tourniquet application and wound packing)
- Chest seal application
- Airway management (NP/OP airways, supraglottic devices, endotracheal intubation)
- Sternal and Tibial intraosseous sites

TCCC MAN[®]

ITEM #	NSN
93-0958	----

Contact us today to learn more or request a **live demonstration.**
www.narescue.com/nar-simulation



INJURIES INCLUDED:

- GSW chest with anterior entry and mid-axillary exit wound
- Deep packable pelvic wound
- Complete below the knee amputation leg (TQ responsive)
- Unstable pelvis with palpable crepitus

SPECIFICATIONS

FEATURE	DETAILS
Power	12V / 5A & 10A NiMH rechargeable batteries (with chargers)
Control Range	Wireless tablet controls up to 12 units
Blood Reservoir	1.2 gallon onboard
Weight	150 lbs.
Durability	Welded steel frame, solid steel joints, rugged silicone skin
Warranty	3-year manufacturer's warranty



NORTH AMERICAN RESCUE[®]
www.NARescue.com • 888.689.6277