

# LASER CUT C-A-T HOLDER

The Laser Cut C-A-T® Holder was specifically designed to allow personnel to place their C-A-T® Tourniquet on their vest or gear for rapid and easy access. It is made of laser cut laminated Cordura nylon. With an innovative construction process this tourniquet holder weighs 50% lighter than traditional nylon TQ holders. The included snap-free PALS/MOLLE-style connectors ensure secure attachment to a duty vest, gear or belt (mounted vertically). To keep your C-A-T® TQ protected from the elements the Laser Cut C-A-T® Holder has an integrated cover with easy-open elastic pull tab.

**NOTE:** Accommodates both Gen 6 and Gen 7 C-A-T® Tourniquets (not included)

Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



## FEATURES/SPECIFICATIONS

- Constructed with laser cut laminated Cordura® nylon
- Easy Pull Elastic Cover
- PALS/MOLLE-style connections
- Weighs less than traditional TQ holders
- Designed & built to easily integrate with new, USDoD laser-cut body armor (PCG3 and MSV) yet "backwards compatible" will easily mount to legacy-type MOLLE vests
- No brass snaps or grommets to wear-out
- Accommodates GEN 6 and GEN 7 C-A-T® tourniquets
- Super durable, yet 50% lighter than traditional nylon TQ holders
- Laser-Cut "TQ" with contrasting background color
- 4 ea. front-mounted, laser cut MAPs (MOLLE Attachment Points) for "stacking". Can easily stack or mount 2 pouches together saving precious plate carrier real-estate
- Deploys easily with one hand, integrated cover with snug, elastic pull-tab
- C-A-T® TQ sold separately

## DIMENSIONS

- Packaged: L 7.25 in. x W 2.75 in. x D 1.5 in.
- Weight: 1.29 oz

## LASER CUT C-A-T HOLDER

ITEM #	COLOR	NSN
30-0297	Black	----
30-0298	Coyote	----
30-5565	Navy	----
30-0556	OD Green	----
30-0299	RGR Green	----
30-0306	Multicam	----
30-0551	Red	----

\* Explore other C-A-T® Holder options at [www.NARescue.com](http://www.NARescue.com)