

TOMMANIKIN LEFT ARM

UNINJURED

The TOMManikin® uninjured left arm features radial pulse, antecubital IV access (with flash), and humeral IO access (procedural only). The simulator's pulse rate, pressure, and regularity controlled via NAR Simulation tablet.

Features 100% silicone construction with metal bone internals. Interchangeable with all TOMManikin® and TAMiKin® simulators.

Learn more about this and other Products with a Mission® at www.NARescue.com



FEATURES/SPECIFICATIONS

- Features radial pulse
- Antecubital IV access
- Humeral IO access (procedural only)
- 100% silicone construction

TOMM LEFT ARM - UNINJURED

ITEM #	NSN
93-0088	----

* Explore more TechLine Trauma products at www.NARescue.com

DIMENSIONS

- L 24 in. x W 12 in. x H 6 in.
- Weight: 10 lbs.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277