

# EZ-IO® G3 EDUCATOR KIT *for Training*

REV113018

## **Complete Educator Kit for the EZ-IO Power Driver...**

EZ-IO® provides the medical professional immediate vascular access to the central circulation within seconds, delivering medications, intravenous fluids and blood products to adult and pediatric patients. With a specially designed cutting IO needle, manual driver, or small power driver, the EZ-IO® allows the clinician complete control — avoiding the use of force. EZ-IO® provides rapid, smooth entry into the bone's medullary cavity, creating an immediate conduit to the central circulation.

The EZ-IO® Intraosseous Infusion System is the leading choice for intraosseous access across multiple health care settings, and is used by 90 percent of US advanced life support ambulances, over half of US Emergency Departments, and the US Military.

Unlike older IO technology, the EZ-IO®'s patented needle tip has four asymmetrical bevels on the stylet, and two cutting tips on the catheter, providing a smooth and accurate insertion. Black indicator lines are located at 5 mm increments from the hub for assessing insertion depth. With EZ-IO®, insertion depth is completely controlled by the clinician.

The G3 Educator Kit comes complete with Training Power Driver, Needle Sets, and Training Bones. All components are packed in a softside Nylon Case for ease of transport and storage.

### **Kit Contents:**

- 1x G3 Training Power Driver
- 1x EZ-IO Adult Training Needle Set
- 1x EZ-IO Pediatric Training Needle Set
- 1x EZ-IO Adult Large Training Needle Set
- 4x Wristband
- 4x EZ-Connect
- 1x G3 Vascular Access Pack
- 4x EZ Stabilizer
- 1x Training Bone-Infant Tibia
- 1x Training Bone-Adult Proximal Humerus
- 1x Training Bone-Adult Tibia

### **EZ-IO® G3 Educator Kit**



Item #	NSN#
30-0074	



***EZ-IO Products are restricted to Federal Purchases only.  
Please contact Customer Service for details.***

*\*Caution: Federal (USA) law restricts this device to sale by or on the order of a practitioner licensed to use or order use of the device or Federal (USA) law prohibits dispensing without a prescription.*



**NORTH AMERICAN RESCUE®**

**www.NARescue.com • 888.689.6277**

**PRODUCTS**  
WITH  **MISSION®**